Description

The HemiCAP® Patello-Femoral Resurfacing Prosthesis incorporates a distal femoral trochlear surface articular component that mates to a fixation stud via a taper interlock, and an all-polyethylene patella component. The prosthesis is intended to be used in cemented arthroplasty.

Materials

Femoral Resurfacing Component:

Cobalt-Chromium Alloy (Co-Cr-Mo)

Undersurface Coating: Titanium (CP Ti)

Fixation Stud: Titanium Alloy (Ti-6Al-4V)

Patella Component:

Ultra-High-Molecular Weight Polyethylene (UHMWPE)

Indications

The HemiCAP® Patello-Femoral Resurfacing Prosthesis is intended to be used in cemented arthroplasty in patients with osteoarthritis limited to the distal patello-femoral joint, patients with a history of patellar dislocation or patellar fracture, and those patients with failed previous surgery (arthroscopy, tibial tubercle elevation, lateral release, etc.) where pain, deformity or dysfunction persists.

Patient selection factors to be considered include:

- 1) Patient's need to obtain pain relief and improve function is significant;
- 2) Patient's tibio-femoral joint is substantially normal;
- 3) Patient exhibits no significant mechanical axis deformity;
- 4) Patient's menisci and cruciates are intact with good joint stability, and good range of motion; and
- 5) Patient's overall well-being is good, including the ability and willingness to follow instructions and comply with activity restrictions.

Contraindications**Absolute contraindications include:**

- 1) Defects that are not localized.
- 2) Inflammatory degenerative joint disease, rheumatoid arthritis, infection, sepsis, or osteomyelitis.
- 3) Patients that have a known sensitivity to materials typically used in orthopedic prosthetic devices or bone cements.

Relative contraindications include:

- 1) Uncooperative patient or patient incapable of following pre-operative and post-operative instructions.
- 2) Metabolic disorders, which may impair the formation or healing of bone; osteoporosis.
- 3) Infections at remote sites, which may spread to the implant site.
- 4) Rapid joint destruction or bone resorption visible on roentgenogram.

- 5) Chronic instability or deficient soft tissues and other support structures.
- 6) Vascular or muscular insufficiency.
- 7) Inadequate skin, musculotendinous or neurovascular system status.

Warnings

Improper selection, placement, positioning, alignment, and fixation of the implant components may reduce the service life of the prosthetic components. Inadequate preparation and cleaning of the implant components mating surfaces may result in improper fixation of the device. Improper handling of the implants can produce scratches, nicks or dents that may have adverse clinical effects on mating joint surfaces. Do not modify implants. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing surgery.

When defining offsets of articular surfaces, care should be taken to ensure that instruments are properly aligned and mated with taper in Fixation Stud. Visually confirm distal tip of Contact Probe is making contact on articular surfaces and free from any soft tissue structures to ensure accuracy. Use light pressure on Contact Probe to slightly indent articular surface at each mapping point, ensuring that the selected implant will be flush or slightly recessed with the articular surface.

Prior to placing implant, carefully trim articular cartilage debris around prepared margin. Remove bone particles and lavage thoroughly. To ensure mechanical interlock of the Fixation Stud and Femoral Resurfacing Component, carefully clean Fixation Stud taper with provided instruments. All drilling or reaming should be done with vigorous lavage to minimize heat effects to adjacent bone and cartilage tissues.

Ensure that care is taken to obtain complete and uniform bone cement coverage at implant site. Unsupported components or unevenly supported components may result in implant failure.

Accepted practices in post operative care should be used. The patient is to be instructed and monitored to ensure a reasonable degree of compliance to post operative instructions and activity restrictions. Excessive activity, impact, and weight gain have been implicated in the reduction of the benefit and service life of prosthetic devices.

Precautions

HemiCAP® Patello-Femoral Resurfacing implants are intended to be fitted and installed with the HemiCAP® Patello-Femoral Resurfacing instrument set. Use of instruments from other systems may result in improper implant selection, fitting, and placement which could result in implant failure or poor clinical outcome. The HemiCAP® instrument set should be regularly inspected for any signs of wear or damage. Do not reuse

implants or disposable instruments. Reuse of single use devices can increase the risk of patient infection and can compromise service life and other performance attributes of the device.

Possible Adverse Effects

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions. Particulate wear debris and mild tissue discoloration from metallic components have been noted in other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
2. Infection or allergic reaction.
3. Loosening, migration or loss of fixation of implant.
4. Fretting and crevice corrosion can occur at the interface between the implant components.
5. Fatigue fracture of the implants as a result of bone resorption around the implant components.
6. Wear and damage to the implant articulating surface.
7. Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.
8. Intraoperative or postoperative bone fracture.
9. Postoperative pain or incomplete resolution of preoperative symptoms.
10. Periarticular calcification or ossification, with or without impediment of joint mobility.
11. Incomplete range of motion due to improper selection or positioning of components.
12. Transient nerve palsy.

Sterility

Metallic prosthetic components are sterilized by exposure to gamma irradiation. Non-metallic prosthetic components are sterilized by gas plasma sterilization. Do not resterilize any components. Do not use components if packaging is opened or damaged. Do not use components if beyond expiration date.

Caution: United States Federal law restricts this device to sale by or on the order of a physician.

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2. Place the **Guide Pin** into a Cannulated Powered Drill and secure at the etch marking on the **Guide Pin**. Advance **Guide Pin** into bone making sure that it is central to the defect. *(It is important to verify that the **Drill Guide** is seated on the curved surface such that all 4 points of contact are established on the articular surface. A normal axis is necessary for proper implant fit.)*

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
7. Clean the taper in the Fixation Stud with Taper Cleaner.
Place Trial Cap into Fixation Stud to confirm correct
depth of Fixation Stud. The height of the Trial Cap must

be flush or slightly below the
existing articular cartilage
surface to avoid the Femoral
Resurfacing Component from
being placed proud or above
the surface of the defect. Adjust
depth if needed using the

Driver to rotate the **Fixation Stud** (rotate clockwise to advance and counterclockwise to retract). Remove **Trial Cap**.

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
8. Place the **Centering Shaft** into taper of **Fixation Stud**. Place **Contact Probe** over **Centering Shaft** and rotate

around shaft. Use light pressure on the **Contact Probe** to ensure proper contact with the articular surface. Read **Contact Probe** to obtain positive (+) superior and inferior offsets, and negative (-) medial and lateral

offsets. Mark each of the identified offsets on the appropriate sizing card. Use the sizing card to record the **maximum superior/inferior** offset and the **minimum medial/lateral** offset.

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
10. Choose the appropriate **Femoral Reamer** based on the **maximum superior/inferior (+)** offset from the sizing card. Confirm selection by matching the color code on the

Femoral Resurfacing Component package with the colored band on the **Femoral Reamer** shaft. Advance **Femoral Reamer** over **Guide Pin** until it contacts the top surface on **Fixation Stud**. (Use lavage

during drilling to prevent possible tissue damage from heat effects). Make sure not to bend the **Guide Pin** during drilling as it may result in **Femoral Resurfacing Component** malalignment.

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13. Prepare and implant the **Patella Component** prior to final placement of the **Femoral Resurfacing Component**.

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
14. Prior to placing the **Femoral Resurfacing Component** on the **Implant Holder** make sure that sufficient suction is

present to hold the device on the distal suction cup. Align the **Femoral Resurfacing Component** on the **Implant Holder**. Make sure to orient the etch marks on the back of the **Femoral Resurfacing**

Component with the etch mark on the handle of the **Implant Holder**. Align the **Femoral Resurfacing Component** with the appropriate offsets. Insert into taper of **Fixation Stud**.

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
1. With knee at 90 degrees flexion, locate the **Alignment Guide** so that the pin fits into the **Fixation Stud**. While observing range of motion, identify target placement of the **Patella Component**

using the pointer on the **Alignment Guide** to transfer **Fixation Stud** central axis. Use slight pressure against the patella so that the pointer on the **Alignment Guide** creates an indentation on distal patella surface.

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
2. Place the **Drill Guide** so that its central axis passes through the indentation on the patella surface. Drill the **Guide Pin** through until it engages the opposite cortex of the patella.

*(It is important to verify that the **Drill Guide** is seated on the curved surface such that all 4*

*points of contact are established on the particular surface.
A normal axis is necessary for proper implant fit).*

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3. Remove the **Drill Guide**. Advance **Circular Scalpel** onto the articular surface to create a cut through the articular surface.

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
4. Place the cannulated **Drill** over **Guide Pin** and drive until the **distal** shoulder of **Drill** is flush to the articular surface. (Use lavage during drilling to prevent possible tissue damage from heat effects). Should the guide pin loosen, use the **Drill** to re-center the **Guide Pin** in the pilot hole and advance into bone.

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
10. Apply a small amount of low-viscosity bone cement onto the underside of the **Patella Component** and quickly place into position.

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
11. Using the **Patella Clamp**, place the *Anatomic or Button* contacting surface against the **Patella Component** and the opposed surface in-line behind the patella. Tighten the **Patella Clamp** until the **Patella Component** is firmly seated in the prepared socket. Leave the **Patella Clamp** in place while the bone cement adequately cures. Remove the **Patella Clamp** and clean out any remaining exposed cement.

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This product is covered by one or more of U.S. Patent Nos. 6,520,964; 6,610,067; 6,679,917; other patents and other patents pending. HemiCAP® is a trademark of ArthroSurface, Inc. U.S.
© 2019 ArthroSurface, Inc. All rights reserved. Printed in the U.S.A.
ArthroSurface, Inc.
MT Promedt Consulting GmbH