Description

The Patellofemoral Wave Arthroplasty Systems incorporate a distal femoral trochlear surface articular component that mates to a taper post via a taper interlock, and an all-polyethylene patella component. The prosthesis is intended to be used in cemented arthroplasty.

Materials

Femoral Articular Component:

Cobalt-Chromium Alloy (Co-Cr-Mo)

Undersurface Coating: Titanium (CP Ti)

Femoral Taper Post:

Titanium Alloy (Ti-6Al-4V)

Patella Component:

Ultra-High-Molecular Weight Polyethylene (UHMWPE)

Indications

Intended to be used in cemented arthroplasty in patients with osteoarthritis limited to the patellofemoral joint, patients with a history of patellar dislocation or patellar fracture, and those patients with failed previous surgery (arthroscopy, tibial tubercle elevation, lateral release, etc.) where pain, deformity or dysfunction persists.

Patient selection factors to be considered include:

- 1) Need to obtain pain relief and improve function is significant.
- 2) Patient's tibio-femoral joint is substantially normal.
- 3) Patient exhibits no significant mechanical axis deformity.
- 4) Patient's menisci and ligaments are intact with good joint stability, and good range of motion.
- 5) Patient age as a potential for early-age-revision of total joint arthroplasty
- 6) Patient's overall well-being, including ability and willingness to follow instructions and comply with activity restrictions.
- 7) Failure of previous conservative treatment options in correcting deformity and achieving pain relief.

Contraindications**Absolute contraindications include:**

- 1) Defects that are not localized.

- 2) Inflammatory degenerative joint disease, rheumatoid arthritis, infection, sepsis, or osteomyelitis.
- 3) Patients that have a known sensitivity to materials typically used in orthopedic prosthetic devices or bone cements.

Relative contraindications include:

- 1) Uncooperative patient or patient incapable of following pre-operative and post-operative instructions.
- 2) Metabolic disorders, which may impair the formation or healing of bone; osteoporosis.
- 3) Infections at remote sites, which may spread to the implant site.
- 4) Rapid joint destruction or bone resorption visible on roentgenogram.
- 5) Chronic instability or deficient soft tissues and other support structures.
- 6) Vascular or muscular insufficiency.
- 7) Inadequate skin, musculotendinous or neurovascular system status.

Warnings

Improper selection, placement, positioning, alignment, and fixation of the implant components may reduce the service life of the prosthetic components. Inadequate preparation and cleaning of the implant components mating surfaces may result in improper fixation of the device. Improper handling of the implants can produce scratches, nicks or dents that may have adverse clinical effects on mating joint surfaces. Do not modify implants. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing surgery.

When defining offsets of articular surfaces, care should be taken to ensure that instruments are properly aligned and mated with taper in Taper Post. Visually confirm distal tip of Contact Probe is making contact on articular surfaces and free from any soft tissue structures to ensure accuracy. Use light pressure on Contact Probe to slightly indent articular surface at each mapping point, ensuring that the selected implant will be flush or slightly recessed with the articular surface.

Prior to placing implant, carefully trim articular cartilage debris around prepared margin. Remove bone particles and lavage thoroughly. To ensure mechanical interlock of the Taper Post and Femoral Component, carefully clean Taper Post taper with

provided instruments. All drilling or reaming should be done with vigorous lavage to minimize heat effects to adjacent bone and cartilage tissues.

Ensure that care is taken to obtain complete and uniform bone cement coverage at implant site. Unsupported components or unevenly supported components may result in implant failure.

Accepted practices in post-operative care should be used. The patient is to be instructed and monitored to ensure a reasonable degree of compliance to post-operative instructions and activity restrictions. Excessive activity, impact, and weight gain have been implicated in the reduction of the benefit and service life of prosthetic devices.

These implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. Their safety in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Precautions

These implants are intended to be fitted and installed with the associated instruments. Use of instruments from other systems may result in improper implant selection, fitting and placement, which could result in implant failure or poor clinical outcome. The instruments should be regularly inspected for any signs of wear or damage.

Surgeon or Physician should discuss general risks and potential complications associated with this and any surgical procedure with the patient prior to patient consent.

Possible Adverse Effects

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions. Particulate wear debris and mild tissue discoloration from metallic components have been noted in other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
2. Infection or allergic reaction.
3. Loosening, migration or loss of fixation of implant.
4. Fretting and crevice corrosion can occur at the interface between the implant components.

5. Fatigue fracture of the implants as a result of bone resorption around the implant components.
6. Wear and damage to the implant articulating surface.
7. Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.
8. Intraoperative or postoperative bone fracture.
9. Postoperative pain or incomplete resolution of preoperative symptoms.
10. Periarticular calcification or ossification, with or without impediment of joint mobility.
11. Incomplete range of motion due to improper selection or positioning of components.
12. Transient nerve palsy.
13. Embolism.

Sterility

Implants and single-use disposable instruments are provided STERILE. Metallic implant components are sterilized by exposure to gamma radiation. Non-metallic implant components (sold separately) are sterilized by gas plasma sterilization. The single-use disposable instruments are sterilized by exposure to gamma radiation. Do not resterilize. Do not use components if packaging is opened or damaged. Do not use components if beyond expiration date. Do not reuse implants or single-use disposable instruments. Reuse of these devices can increase the risk of patient infection and can compromise service life and other performance attributes of the device(s).

Caution

United States Federal law restricts this device to sale by or on the order of a physician.

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4. Select the **Guide Block** that corresponds with the S/I curvature value and place onto the trochlear groove. Insert the central portion of the **Guide Block** into the reamed area and secure the **Guide Block** onto the femur using the **2.0mm Guide Pins**.

In knees with a dysplastic or flattened trochlea, the **Guide Block** may not sit flush to the reamed area. Use an osteotome or burr to create slots in the bone to accept the proximal and distal feet of the **Guide Block**. Creating these slots will allow the **Guide Block** to sit flush to the reamed area.

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5. Assemble the **Outer Reamer** into the **Guide Bushing**. Secure the **Guide Bushing** into the inferior **Guide Block** bore. Advance the **Outer Reamer** into the bone until the depth mark on the reamer shaft is reached. Remove assembly and repeat reaming through the superior **Guide Block** bore. It is critical to keep the **Guide Block** stable during reaming. Repeat the same sequence for the cannulated **Kahuna Tooth Cutter**.

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8. Assemble the **Trial Handle** onto the **Sizing Trial** that corresponds with the S/I curvature value. Place the **Sizing Trial** into the prepared site. Confirm the fit of the **Sizing Trial** so that all margins are congruent or recessed to the edge of the surrounding articular surface. Trim the transition areas between the reamed surfaces to ensure the **Sizing Trial** is fully seated. Fix the **Sizing Trial** in place using the **2.0mm Guide Pins**.

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
12. Apply a small amount of low-viscosity bone cement into the **Taper Post** bone tunnel. Place the **Taper Post** into the Morse taper of the **Trial Handle** and attach to the **Sizing Trial**. Place the **Hex Driver** through the **Trial Handle** and advance the **Taper Post** until the stop on the shaft of the **Hex Driver** comes in contact with the back of the **Trial Handle**. Place the **Placement Gauge** into the **Sizing Trial** to ensure that the **Taper Post** is at proper depth to engage the **WaveKahuna™ Femoral Component**.

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13. An alternative approach to fixation is to pre-assemble the **Threadless Stud** to the **WaveKahuna™ Femoral Component**. Be sure to protect the articular face of the **Kahuna™ Femoral Component** by using slight impaction with the mallet to seat the morse taper of the **Threadless Stud** onto the **Kahuna™ Femoral Component**.

NOTE: Prepare the **Patella Component** implant site prior to the final placement of the **WaveKahuna™ Femoral Component**.

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CE
0459

Rx ONLY

This product is covered by one or more of U.S. Patent
Nos. 6,520,964; 6,610,067; 6,679,917; other patents and
other patents pending. HemiCAP® is a trademark of
Arthrosurface, Inc. U.S.
© 2022 Arthrosurface, Inc. All rights reserved.
Printed in the U.S.A.

PN 2001-2120 REV D

Patellofemoral Wave Arthroplasty Systems
Instructions for Use

*PF Wave® and WaveKahuna® are marketed and sold
under the trademarks PFXL™ and PF Kahuna™ in
the EU and UK*

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MT Promedt Consulting GmbH