Parcus V-LoX™ Titanium Suture Anchor

Important Product Information

Produktinformationsblatt

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Parcus V-LoX™ Titanium Suture Anchor (English)

1. Indications

The Parcus V-LoX™ Titanium Suture Anchors are indicated for attachment of soft tissue to bone. This product is intended for the following indications:

Shoulder Rotator Cuff Repair, Acromioclavicular Separation Repair, Bankart Lesion Repair, Biceps Tenodesis, Capsular Shift or Capsulolabral Reconstruction, Deltoid Repair, SLAP Lesion Repair.

Knee Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Extra Capsular Reconstruction, Iliotibial Band Tenodesis, Patellar Ligament and Tendon Avulsion Repair.

Foot/Ankle Lateral Stabilization, Medial Stabilization, Midfoot Reconstruction, Achilles Tendon Repair, Hallux Valgus Reconstruction, Metatarsal Ligament Repair.

Elbow Tennis Elbow Repair, Biceps Tendon Reattachment.

Hand/Wrist Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction, TFCC.

Hip Acetabular Labral Repair

2. Contraindications

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the instructions for use or using the device for an indication other than those listed.

3. Adverse Effects

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings

- A. This product should only be used by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Anchor size selection should be made with care taking into consideration the quality of the bone in which the anchor is to be placed. Osteopenic bone poses fixation challenges which may be addressed by larger diameter V-LoX Titanium Suture Anchors.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- E. Polyethylene polyblend sutures, as used in this device, elicit a minimal inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Polyethylene polyblend suture is not absorbed, nor is any significant change in tensile strength retention known to occur in vivo.
- F. The patient should be advised of the use and limitations of this device.
- G. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- H. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- I. This device must never be re-sterilized.
- J. Appropriate instrumentation should be used to implant this device.
- K. This device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment. Use of MR technology in the presence of devices of this nature may cause magnetically induced displacement forces and torques, radio frequency heating and image artifacts. Standard MRI screening guidelines for post-operative patients should be followed.

5. **Packaging and Labeling**

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. **Material Specifications**

The V-LoX Titanium Suture Anchor is supplied with or without needles attached to the high-strength, braided, polyethylene polyblend sutures, The anchor material is Ti-6Al-4V ELI (ASTM F136).

7. **Sterilization**

The V-LoX Titanium Suture Anchor and driver are supplied in sterile packaging. The contents are sterilized by EO gas. These products must never be re-sterilized.

8. **Storage**

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. **Instructions for Use**

- A. Identify bone of sufficient quantity and quality into which the anchor is to be placed.

Note: In cases where bone quality is suspect and the 4.5mm V-LoX Titanium Suture Anchor does not provide the desired fixation, the larger diameter 5.5mm or 6.5mm V-LoX Titanium Suture Anchor may provide an alternative solution.

- B. Select the desired configuration (with or without needles) and aseptically open the package containing the V-LoX Titanium Suture Anchor.
- C. Taking care to create an entry to the desired placement site as close to perpendicular as possible, place the V-LoX Titanium Suture Anchor and driver assembly through a suitable arthroscopic cannula, percutaneously, or in an open case, directly into the wound. Using a suitable mallet, gently strike the driver handle until the threads of the anchor contact the bone.
- D. Turn the driver handle clockwise until the circumferential laser etched line on the distal end of the driver shaft is flush with the surrounding bone. If present the longitudinal lines on the driver shaft are used to indicate the orientation of the suture passing around the post within the anchor. Placing the longitudinal lines in a plane perpendicular to the targeted tissue offers the best orientation to allow the suture to slide through the anchor.
- E. If the V-LoX Titanium Suture Anchor with Needles is selected, release the suture loops from the back of the handle and disengage driver from anchor. The needles will now be loose and may be removed.
- F. Two or three strands of polyethylene polyblend high strength suture either with or without needles are attached to the anchor. Note that each strand is different color in order to aid in identification.
- G. Pass the sutures through the targeted tissue in the user preferred arthroscopic manner. In open cases, use the needles to pass the suture.
- H. Complete the tissue repair by tying knots of the surgeon preferred configuration and cutting the suture tails above the knots.

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