Anika PEEK CF Interference Screw Driver
Destornillador de interferencia PEEK CF Anika
Tournevis pour vis d'interférence Anika PEEK CF
Anika PEEK CF Interferenzschrauber
Cacciavite ad interferenza Anika PEEK CF
Chave de Parafuso de Interferência PEEK CF Anika
Κατσάβιδι Παρεμβολής Anika PEEK CF
„Anika PEEK CF Interference“ sraigčių suktuvas
Anika PEEK CF interferensskrutrekker
Anika PEEK CF skruvmejsel för interferensskruvar
Șurubelniță pentru șurub de interferență PEEK CF Anika
Wkrętał interferencyjny z polieteroeteroketonu i włókien węglowych (PEEK CF) Anika
Отвертка за поддържащи винтове Anika PEEK CF

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30350 Rev 7
EIT: 07/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 ANIKA
www.anika.com
(941) 755-7965

 2797
EC REP
MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

Anika PEEK CF Interference Screw Driver (EN)

1. Intended Use:

The Anika PEEK CF Interference Screw Driver is designed to insert and remove Anika PEEK CF Interference Screws. The shafts feature a 1.5mm cannulation and are designed to fit into the Anika Ratcheting Handle and may fit into other handles that utilize a Hudson chuck.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. This product should not be used to implant any device other than Anika PEEK CF Interference Screws.
- C. Do not use this product if it is bent or there is visible damage to the shaft or drive section.

3. Material:

This device is manufactured from a Stainless-Steel Alloy. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

6. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹Drying times vary according to load size and should be increased for larger loads.

²This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

7. Directions for use:

- A. Connect the Anika PEEK CF Interference Screw Driver to the Anika Ratcheting Handle by pulling the locking collar of the handle towards its proximal end and inserting the Anika PEEK CF Interference Screw Driver until it is fully seated. It will be necessary to align the flats of the driver with the slots of the handle to allow the locking collar to slide distally, securing the driver to the handle. (See Anika Ratcheting Handle Directions for Use for details on attaching accessories).
- B. Place an Anika PEEK CF Interference Screw onto the drive portion of the driver until it is fully engaged.
- C. Place the screw and driver over an appropriately placed Anika 1.5mm Nitinol Guide Wire to guide the screw into the insertion site.
- D. Rotate the driver clockwise until the screw has reached the desired depth. A counterclockwise rotation of the driver will remove the screw.
- E. Remove the driver from the screw followed by removal of the guide wire.

8. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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