Anika 3.5mm Hex Driver with Handle
Destornillador hexagonal de 3,5 mm Anika con mango
Tournevis hexagonal avec poignée 3,5mm Anika
Anika 3,5mm-Sechskantschrauber mit Griff
Cacciavite esagonale Anika da 3,5 mm con impugnatura
Chave Sextavada de 3,5 mm com Alça Anika
Εξάγωνο Κατσαβίδι 3.5mm με λαβή Anika
„Anika“ 3,5 mm šešiabriaunis suktuvus su rankena
Anika 3,5 mm sekskanttrekker med håndtak
Anika 3,5 mm sexkantsskruvmejsel med handtag
Șurubelniță hexagonală Anika de 3,5 mm cu mâner
Wkrętak sześciokątny o rozmiarze 3,5 mm z uchwytem Anika
Шестограм с дръжка 3,5 мм Anika

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30385 Rev 8
Eff. 07/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 ANIKA  2797
www.anika.com
(941)755-7965

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Anika 3.5mm Hex Driver with Handle (EN)

1. Intended Use:

The Anika 3.5mm Hex Driver with Handle is designed to insert and remove Anika Titanium Interference Screws. The driver features a 1.5mm cannulation and is compatible with Anika 1.5mm Guide Wires.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. This product should not be used to implant any device other than Anika Titanium Interference Screws.
- C. Do not use this product if it is bent or there is visible damage to the shaft or drive section.

3. Material:

This device is comprised of a Stainless-Steel Alloy shaft and an Aluminum handle. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

6. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

7. Directions for use:

- A. Place an Anika Titanium Interference Screw onto the drive portion of the driver until it is fully engaged.
- B. Place the screw and driver over an appropriately placed Anika 1.5mm Guide Wire to guide the screw into the insertion site.
- C. Rotate the driver clockwise until the screw has reached the desired depth. A counterclockwise rotation of the driver will remove the screw.

Remove the driver from the screw, followed by removal of the guide wire.

8. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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