Anika 5.5 – 6.5mm Punch/Tap with Handle
Anika 5,5 – 6,5 mm Punzón/terrajá con mango
Mucotome/taraud 5.5 – 6.5mm avec poignée Anika
Anika 5,5 – 6,5mm-Stanzinstrument/Gewindebohrer mit Griff
Punzone Anika da 5,5– 6,5 mm con impugnatura
Puncionador/Percussionador com Alça de 5,5 – 6,5 mm Anika
Διατρητής/Βύσμα με Λαβή 5,5 – 6,5mm Anika
„Anika“ 5,5–6,5 mm perforatorius (sriegiklis)
Anika 5,5–6,5 mm dor/gjengetapp med håndtak
Anika 5,5–6,5mm stans/tapp med handtag
Dorn/burghiu Anika de 5,5-6,5 mm cu mâner
Wycinak/Gwintownik o rozmiarze 5,5 – 6,5 mm z uchwytem Anika
Перфоратор/накрайник с дръжка 5,5 - 6,5 мм Anika

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30387 Rev 7
EUF. 07/24

Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 **ANIKA**  2797
www.anika.com
(941)755-7985

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Anika 5.5 – 6.5mm Punch/Tap with Handle (EN)

1. Intended Use:

The Anika 5.5 – 6.5mm Punch/Tap with Handle is designed to prepare bone for the implantation of the Anika 5.5mm and 6.5mm V-Lox™ PEEK CF and Twist PEEK Suture Anchors.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. This product should not be used to implant any device other than the Anika V-Lox PEEK CF or Twist PEEK Suture Anchors.
- C. Do not use this product if it is bent or there is visible damage to the tip or thread form.

3. Material:

This device is manufactured from a Stainless-Steel Alloy shaft and an Aluminum handle. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

5. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10-6 sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

6. Directions for use:

- A. Identify the site for anchor insertion and gently tap the Anika Punch/Tap into the bone until the threads of the Punch/Tap can be engaged.
- B. Rotate the Punch/Tap clockwise into the bone until the top of the laser line is flush with the bone surface.
- C. Carefully remove the Punch/Tap by rotating the handle in a counter-clockwise direction until its threaded portion is clear of the bone.
- D. Insert an Anika V-Lox PEEK CF or Twist PEEK Suture Anchor into the prepared hole.

7. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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