Anika Cannulated Drills and Headed Reamers
Fresas canuladas y escariadores con cabeza Anika
Forets canulés et alésoirs à tête Anika
Punte cannulate e alesatori con testa Anika
Anika kanülierte Bohrer und Kopfreibahlen
Brocas Canuladas e Escareadores com Cabeça Anika

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Informazioni importanti sul prodotto
Wichtige Produktinformationen
Informações Importantes sobre o Produto

Instructions for Use
Instrucciones de uso
Mode d'emploi
Istruzioni per l'uso
Gebrauchsanweisung
Instruções de uso

30470 Rev 9
Eff. 08/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 **ANIKA**
www.anika.com
(941)755-7965

 2797

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Anika Cannulated Drills and Headed Reamers (EN)

1. Intended Use

Anika Cannulated Drills and Headed Reamers are used to create tunnels and sockets in bone. This product is designed to be used in conjunction with a 2.4mm guide pin.

2. Warnings:

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. Caution: Common risks consistent with surgery may include: bone/tissue damage, revision surgery, soft tissue irritation, infection, allergic reactions.
- C. Never place this product over a bent or damaged guide pin.
- D. Placing off-axis loads on the guide pin by deflecting the drill/reamer may cause damage to both the drill/reamer and the guide pin.
- E. Never use a dull or damaged drill/reamer.
- F. This product is not intended to be left in the body.
- G. Be sure to remove any debris from the cannulation of the drill/reamer during the cleaning process.

3. Material:

These products are manufactured from a stainless-steel alloy. The material used in the manufacture of the devices that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

Instruments intended to be reused should be inspected prior to each use to assess damage and wear. Signs of damage and wear include surface discoloration, surface corrosion, surface damage, and general bending, cracking, and fractures. Cutting, drilling, and other instruments intended to have sharp edges will become dull over time, diminishing their ability to function as designed. Review sharp instruments for dullness or damage at the cutting/sharp feature. If the instrument does not pass inspection, do not use the instrument, and notify an Anika representative.

5. Storage:

Products must be stored in a dry place.

6. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

7. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10^{-6} sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

After sterilization aseptically transfer device into the sterile field in accordance with the institution's policy.

8. Instructions for use:

- Aseptically transfer device into sterile field and use a purpose designed drill motor or wire driver to place a guide pin in the desired position.
- Securely attach the drill/reamer into a drill chuck.
- Place the drill/reamer over the guide pin.
- Make sure that the drill motor is set to rotate the drill/reamer in a clockwise direction.
- Prudently advance the drill/reamer over the guide pin, ensuring to maintain a coaxial orientation between the pin and the drill/reamer.
- After the desired depth has been reached, reverse the drill motor and remove the drill/reamer from the bone.

9. Device Disposal:

Devices that have evidence of damage must be properly disposed of in accordance with the institution's policy.

10. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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