Parcus Femoral Offset Guide

Important Product Information

Produktinformationsblatt
Información Importante sobre el Producto
Informations Importantes sur le Produit
Informazioni Importanti sul Prodotto
Informações Importantes sobre o Produto

Directions for Use

Gebrauchsanleitung
Instrucciones de uso
Mode d'emploi
Istruzioni per l'uso
Instruções para o Uso

30472 Rev 6
Eff. 09/22

 Mfd by:
Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243

ANIKA

www.anika.com
(941)755-7965

2797

EC REP

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Parcus Femoral Offset Guide

1. Indications:

The Parcus Femoral Offset Guides are used to assist in placing a 2.4mm guide pin for creating the femoral tunnel or socket in ACL reconstructions.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. No component of this instrument is intended to be left in the body.
- C. Never use this device if it has been bent or damaged in any way.
- D. Extending the knee with this instrument in position may result in damage to the instrument.
- E. Take care to remove all debris from the cannulation of the instrument after use.

3. Material:

This product is manufactured from stainless steel and aluminum alloys. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to the cannulation and other areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

6. Sterilization:

This product is a non-sterile item that must be cleaned and sterilized prior to use. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 minutes	30 minutes ¹
Pre-vacuum	132°C / 270°F	4 minutes	

¹ Drying times vary according to load size and should be increased for larger loads

7. Directions for use:

- A. Select the proper offset. The Femoral Offset Guides are sized in various millimeter offsets that place a guide pin anterior to the tang of the guide. For example, a 7mm aimer will leave a 2mm wall of bone if a 10mm reamer is used.
- B. Insert the aimer into the knee through a tibial tunnel or arthroscopic portal. It may be necessary to flex the knee to place the guide in the desired position.
- C. Place the tang of the Femoral Offset Guide at the “over the top” position on the femur.
- D. Advance a 2.4mm guide pin through the cannulation in the guide into or through the bone.
- E. Remove the Femoral Offset Guide and use a cannulated drill or reamer over the guide pin to create the tunnel or socket in the targeted bone.

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