Anika Graft Station
Estación de injerto Anika
Système de préparation de greffe Anika
Anika Transplantatstation
Stazione d'innesto Anika
Estação de Enxerto Anika
Σταθμός Μοσχεύματος Anika
„Anika“ transplantantų stotis
Anika transplantasjonsstasjon
Anika transplantatstation
Stację de preparare grefă Anika
Stacja do przeszczepów Anika
Станция за присадки Anika

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30474 Rev 8
Eff. 08/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

www.anika.com
(941)755-7965

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Anika Graft Station (EN)

1. Intended Use:

The Anika Graft Station is used to assist in the preparation of graft tissue and material used in the repair of ligaments and tendons.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. No component of this instrument is intended to be left in the body.

3. Material:

This product is manufactured from stainless steel, aluminum alloys, and autoclavable plastics.

4. Inspection:

Instruments intended to be reused should be inspected prior to each use to assess damage and wear. Signs of damage and wear include surface discoloration, surface corrosion, surface damage, and general bending, cracking, and fractures. Cutting, drilling, and other instruments intended to have sharp edges will become dull over time, diminishing their ability to function as designed. Review sharp instruments for dullness or damage at the cutting/sharp feature. If the instrument does not pass inspection, do not use the instrument, and notify an Anika representative.

5. Storage:

Products must be stored in a dry place.

6. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

7. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10^{-6} sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

After sterilization aseptically transfer device into the sterile field in accordance with the institution's policy.

8. Directions for use:

- A. Aseptically transfer device into sterile field and assemble the Anika Graft Station by placing the Graft Grasper and Tensioner into the keyhole slot in the base. Orient the Graft Grasper and Tensioner so that the large tensioning knob is furthest away from the graft post. Select the Fixed Clamp or GFS (Graft Fixation System) Post and place it in one of the two slots opposite the Graft Grasper and Tensioner and lock it to the base by turning the thumb screw at the base of the post.
- B. Adjust the Graft Grasper and Tensioner so that it can capture the graft by moving it along the slot in the Graft Station Base. Lock the Graft Grasper and Tensioner in position by turning the large knob nearest the base, in a clockwise direction.
- C. The graft can be secured to the Graft Station by looping it mid-substance around the GFS Post and securing the ends of the graft by using the clamps on the Graft Grasper and Tensioner.
- D. Once the graft has been prepared, the large knob on the Graft Grasper and Tensioner can be rotated in a clockwise manner to place tension on the graft. The calibration marks on the Graft Grasper and Tensioner can be used to determine the tension placed on the graft. The calibration marks represent the tension in pounds. Conversion to kilograms is as follows:

10 lbs = 4.5 kg

15 lbs = 6.8 kg

20 lbs = 9.1 kg

9. Device Disposal:

Devices that have evidence of damage must be properly disposed of in accordance with the institution's policy.

10. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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