Anika Threaded Dilator
Dilatador Roscado Anika
Dilatateur fileté Anika
Anika Gewindedilator
Dilatatore filettato Anika
Dilatador Rosqueado Anika
Διαστολέας με Σπείρωμα Anika
„Anika“ srieginis plėstuvė
Anika gjenget dilatator
Anika gängad dilatator
Dilatator filetat Anika
Rozszerzacz gwintowany Anika
Дилататор с резьба Anika

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Instructions for Use
Instrucciones de uso
Mode d'emploi
Gebrauchsanweisung
Istruzioni per l'uso
Instruções de uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Användarinstruktioner
Instrucțiuni de utilizare
Instrukcja użytkowania
Инструкции за употреба

30635 Rev 7
Eff. 08/24

 **Parcus Medical, LLC**
6423 Parkland Dr.
Sarasota, FL 34243 USA

 **ANIKA**
www.anika.com
(941) 755-7965

CE 2797
EC REP
MDSS GmbH
Schiffgraben 41
30175 Hanover, Germany

Anika Threaded Dilator (EN)

1. Intended Use:

The Anika Threaded Dilator is designed to create a starter hole to ease the insertion of Interference Screws in Bone-Tendon/Ligament-Bone ACL Reconstructions. The Threaded Dilator is cannulated and is compatible with 1.5mm or less Guide Wires.

2. Warnings:

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. Do not use this product if it is bent or there is visible damage to the shaft or drive section.
- C. This device is only to be used when sufficient bone quality and implant/graft size exists such that fixation is not compromised.
- D. This device is not intended to be used as an implant.
- E. Placing off-axis loads on the threaded dilator may cause damage to the device.
- F. Placing off-axis loads on the guide pin by deflecting the threaded dilator may cause damage to both the threaded dilator and the guide pin.
- G. Caution: Common risks consistent with surgery may include: bone/tissue damage, revision surgery, soft tissue irritation, infection, allergic reactions.

3. Material:

This device is manufactured from a Stainless-Steel Alloy and may include an aluminum alloy handle. The materials used in the manufacture of this device are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

Instruments intended to be reused should be inspected prior to each use to assess damage and wear. Signs of damage and wear include surface discoloration, surface corrosion, surface damage, and general bending, cracking, and fractures. Cutting, drilling, and other instruments intended to have sharp edges will become dull over time, diminishing their ability to function as designed. Review sharp instruments for dullness or damage at the cutting/sharp feature. If the instrument does not pass inspection, do not use the instrument, and notify an Anika representative.

5. Storage:

Products must be stored in a dry place.

6. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

7. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum

sterilization parameters that have been validated by Anika Therapeutics to provide a 10^{-6} sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

After sterilization aseptically transfer device into the sterile field in accordance with the institution's policy.

8. Instructions for use:

- Aseptically transfer device into sterile field and place the Threaded Dilator over an appropriately placed 1.5mm Guide Wire to guide the dilator into the insertion site.
- Rotate the Threaded Dilator clockwise until the dilator has reached the desired depth.
- Rotate the Threaded Dilator counterclockwise to remove the dilator.
- Place the desired Interference Screw over the Guide Wire and insert per applicable Instructions for Use.

9. Device Disposal:

Devices that have evidence of damage must be properly disposed of in accordance with the institution's policy.

10. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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