Parcus Braid™
Non-Absorbable Surgical Suture, Infinity Loop,
and Parcus Suture Tape

Important Product Information

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1. Indications:

Parcus Braid products are indicated for use in general soft tissue approximation and/or ligation and the use of allograft tissues for orthopedic surgeries. Parcus Braid products are intended for one-time use only and are not to be re-sterilized.

2. Contraindications:

None known.

3. Adverse Effects:

- A. Adverse effects associated with the use of this device include: wound dehiscence, calculi formation when prolonged contact with salt solutions occur, enhanced bacterial infectivity, minimal acute inflammatory reaction, and transitory local irritation.
- B. Parcus Braid products elicit a minimal acute inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Parcus Braid products are not absorbed, nor is there any significant change in tensile strength retention known to occur in vivo.

4. Warnings:

- A. CAUTION: Federal (USA) law restricts this device to sale and use by, or on the order of, a physician.
- B. As with any foreign body, prolonged contact of any suture products with salt solutions may result in calculus formation.
- C. Do not resterilize. Discard opened or unused sutures.
- D. Users should be familiar with surgical procedures and techniques involving suture products before employing Parcus Braid products are used for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

5. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.
- C. The sutures are supplied sterile, in pre-cut lengths, non-needed or affixed to stainless steel needles. Parcus Braid products are available in quantities of one, six or twelve per box.

6. Material Specifications:

Parcus Braid Sutures, Infinity Loop and Parcus Suture Tape are non-absorbable, sterile, surgical suture products composed of ultra high molecular weight polyethylene (UHMWPE). The multifilament braided sutures are available uncoated. Parcus Braid sutures are inert and elicit only minimal local tissue reaction. The suture braid is available undyed (white), dyed blue, dyed black or with trace filaments of black nylon, blue PET, or green PET. Parcus Braid products may be provided with stainless steel needles or tipped ends using cyanoacrylate. Parcus Braid Suture is available in multiple color configurations ranging in sizes from USP 6-0 through USP 7. The Parcus Braid Infinity Loop is available in sizes from USP 2 through USP 5. Parcus Braid Suture Tape is available in widths from 1mm through 5mm. Parcus Braid products meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopeia (EP) for non-absorbable surgical sutures except for diameter and where indicated on the individual unit package. These devices are MR Safe.

7. Sterilization:

Parcus Braid Sutures, Infinity Loop and Parcus Suture Tape are supplied in sterile packaging. The contents are sterilized by EO gas. These products must never be re-sterilized.


8. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Instructions for Use:

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon.

Definition of symbols:

Symbol	Definition
	Manufacturer
	Authorized representative in the European Community
	This product meets the essential requirements of the Medical Device Directive 93/42 EEC as reviewed by NB 2797.
	Use-by-date
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Do not re-sterilize
	Do not use if package is damaged
	Do not re-use
	Consult instructions for use
	Consult electronic instructions for use
	MR safe
	Medical device

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