Anika Reusable Knotless PEEK CF Suture Anchor Driver
Destornillador de anclaje de sutura PEEK CF reutilizable sin nudos Anika
Tournevis à usage multiple pour Ancre de suture sans nœud PEEK CF Anika
Anika Wiederverwendbarer Schrauber für Verankerungen von knotenfreien PEEK CF Fäden
Cacciavite per ancora di sutura Anika PEEK CF senza nodi riutilizzabile
Chave de Âncora de Sutura Reutilizável Knotless PEEK CF Anika
Επαναχρησιμοποιήσιμος Οδηγός Αγκιστρών Ραμμάτων Χωρίς Κόμπους Anika PEEK CF
„Anika“ daugkartinio naudojimo bemazgių PEEK CF siūlų inkarų suktuvas
Anika gjenbrukbar Knotless PEEK CF suturankertrekker
Anika återanvändbara knutfria PEEK CF drivanordning för suturankare
Șurubelniță Anika pentru ancore de sutură PEEK CF fără noduri, reutilizabile
Wkrętak wielokrotnego użytku do bezwęzłowej kotwicy do szwów z PEEK CF Anika
Водач за многократна употреба за шевни анкери без възли PEEK CF Anika

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30892 Rev 6
Eff. 08/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 ANIKA
www.anika.com
(941)755-7965

 2797
EC REP
MDSS GmbH
Schiffgraben 41
30175 Hanover, Germany

Anika Reusable Knotless PEEK CF Suture Anchor Driver (EN)

1. Intended Use:

Anika Knotless PEEK CF Suture Anchor Drivers are provided unsterile as accessories to Knotless PEEK CF Suture Anchors. They are used to assist in the insertion of Knotless PEEK CF Suture Anchors in bone.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. Never use a bent driver.
- C. This product is not intended for use as an implant.
- D. Loading the driver in an off-axis manner may result in breakage of the implant or bending of the driver.

3. Material:

This device is manufactured from stainless steel and aluminum alloys. The materials used in the manufacture of this device are radiopaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected consult the manufacturer for guidance.

5. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

6. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

7. Directions for use:

- A. Aseptically open the package containing the Knotless PEEK CF Suture Anchor.

- B. Assemble the suture anchor to the appropriately sized driver by inserting the tip of the driver into the socket of the anchor and pressing until the driver is firmly attached to the anchor.

Warning: Off-axis or insufficient seating of the anchor to the driver may result in breakage of the implant.

- C. Pass sutures and insert implant into bone using the preferred user technique as described in implant DFU.
D. Once implant is fully inserted, remove driver by pulling firmly and twisting slightly until the driver is free from the anchor.

8. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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