Anika Switching Stick
Barra de conmutación Anika
Switching Stick Annika
Anika Instrumentenwechsler
Asta di sezionamento Anika
Bastão de Comutação Anika
Κοντάρι Εναλλαγής Anika
„Anika“ keitimo lazdelė
Anika skiftepinne
Anika bytesstav
Tijă de interschimbare Anika
Szyft przełączający Anika
Стик за превключване Anika

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30896 Rev 6
Eff. 08/24

Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

ANIKA
www.anika.com
(941)755-7965

2797
EC REP
MDSS GmbH
Schiffgraben 41
30175 Hanover, Germany

Anika Switching Stick (EN)

1. Intended Use:

The Anika Switching Stick is designed to assist in the movement of cannulated instruments in and out of the various joints and to maintain the orientation and location of holes prepared in bone. The Switching Stick is designed to work with various Drill Guides and cannulas commonly used in arthroscopic procedures.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. Do not use this product if it is bent or there is visible damage.
- C. This product is not intended for use as an implant.

3. Material:

This device is manufactured from a Stainless-Steel Alloy. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

6. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination.

7. Directions for use:

- A. Place the Switching Stick through a skin portal and pass an instrument, with a cannulation large enough to accommodate its diameter, over the Switching Stick and into the joint.

- B. Conversely, the Switching Stick may be used to exchange a previously placed cannulated instrument, with another, by placing the Switching Stick through the first instrument. The first instrument is extracted, leaving the Switching Stick behind. The second instrument is then passed over the Switching Stick and into the joint.
- C. The Switching Stick may be used as a marker to assist the surgeon in marking the location and directional orientation of holes drilled in bone. After using a drill bit through a guide to create a hole in the targeted bone, the drill is removed from the guide and replaced with the Switching Stick. The user should ensure that the Switching Stick is inserted into the full depth of the hole. The guide is then removed leaving the Switching Stick behind to preserve the integrity of the hole and mark its location.

8. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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