Parcus Extended Shaft MiTi Suture Anchor

Important Product Information

Produktinformationsblatt
Información Importante sobre el Producto
Informations Importantes sur le Produit
Informazioni Importanti sul Prodotto
Informações Importantes sobre o Produto

Directions for Use

Gebrauchsanleitung
Instrucciones de uso
Mode d'emploi
Istruzioni per l'uso
Instruções para o Uso

30962 Rev 5
Eff. 08/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 **ANIKA**  2797
www.anika.com
(941) 755-7985
MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

Parcus Extended Shaft MiTi Suture Anchor

1. Indications:

The Parcus MiTi Suture Anchors are indicated for attachment of soft tissue to bone. This product is intended for the following indications:

<u>Shoulder</u>	Rotator Cuff Repair, Acromioclavicular Separation Repair, Bankart Lesion Repair, Biceps Tenodesis, Capsular Shift or Capsulolabral Reconstruction, Deltoid Repair, SLAP Lesion Repair.
<u>Knee</u>	Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Extra Capsular Reconstruction, Iliotibial Band Tenodesis, Patellar Ligament and Tendon Avulsion Repair.
<u>Foot/Ankle</u>	Lateral Stabilization, Medial Stabilization, Midfoot Reconstruction, Achilles Tendon Repair, Hallux Valgus Reconstruction, Metatarsal Ligament Repair.
<u>Elbow</u>	Tennis Elbow Repair, Biceps Tendon Reattachment.
<u>Hand/Wrist</u>	Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction, TFCC.

2. Contraindications:

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected, should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the directions for use or using the device for an indication other than those listed.

3. Adverse Effects:

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings:

- A. The Parcus MiTi Suture Anchors have not been evaluated for safety and compatibility in the MR environment. The Parcus MiTi Suture Anchors have not been tested for heating or migration in the MR environment.
- B. This product should only be used by or on the order of a physician.
- C. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- D. Anchor size selection should be made with care taking into consideration the quality of the bone into which the anchor is to be placed. Osteopenic bone poses fixation challenges which may be addressed by larger diameter suture anchors as the size of the anchor placement site dictates.
- E. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- F. Polyethylene polyblend sutures, as used in this device, elicit a minimal inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Polyethylene polyblend suture is not absorbed, nor is any significant change in tensile strength retention known to occur in vivo.
- G. The patient should be advised of the use and limitations of this device.

- H. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity with the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
 - I. This device must never be reused. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
 - J. This device must never be re-sterilized.
 - K. Appropriate instrumentation should be used to implant this device.
 - L. This device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment. Use of MR technology in the presence of devices of this nature may cause magnetically induced displacement forces and torques, radio frequency heating and image artifacts. Standard MRI screening guidelines for post-operative patients should be followed.
- 5. Packaging and Labeling:**
- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
 - B. Please contact Parcus Medical Customer Service to report any package damage or alterations.
- 6. Material Specifications:**
- The MiTi Suture Anchor is supplied with high-strength, braided, polyethylene polyblend sutures. Contrasting suture colors are used in suture anchors supplied with multiple sutures. The anchor material is Ti-6Al-4V ELI (ASTM F136). The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.
- 7. Sterilization:**
- The MiTi Suture Anchor and driver assembly is supplied sterile. These products must never be re-sterilized.
- 8. Storage:**
- Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.
- 9. Directions for Use**
- A. Identify bone of sufficient quantity and quality into which the anchor is to be placed.

Note: In cases where bone quality is suspect and a MiTi Suture Anchor does not provide the desired fixation, a larger diameter suture anchor may provide an alternative solution.

- B. Arthroscopic cannulae, are typically used in conjunction with these devices. Always try to approach the targeted site as close to perpendicular as possible. When using the Drill Guide, use the obturator supplied with the Drill Guide to aid in its passage through tissue. Place the distal tip of the guide on the targeted bone and insert the Suture Anchor through its lumen.
- C. Place the MiTi Suture Anchor, as described above, and turn the driver handle clockwise until the circumferential laser etched line on the distal end of the driver shaft reaches the end of the drill guide window. The orientation of the sutures passing through the implant may be visualized as they are exterior to the shaft. Placing the sutures in a plane perpendicular to the targeted tissue offers the best orientation to allow the suture to slide through the anchor.
- D. One or two strands of polyethylene polyblend high strength suture are attached to the anchor.
- E. Pass the sutures through the targeted tissue in the user preferred manner.
- F. Complete the tissue repair by tying knots of the surgeon preferred configuration and cutting the suture tails above the knots.

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