**Directions for Use
Parcus Chondral Picks**

1. Indications:

The Parcus Chondral Picks are designed to create a socket/defect in bone.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. Do not use this product if it is bent or there is visible damage to the shaft or drive section.
- C. This device is not intended to be used as an implant.

3. Material:

This device is manufactured from a Stainless Steel Alloy and may include an aluminum alloy handle. The materials used in the manufacture of this device are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

6. Sterilization:

This device is not provided sterile. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a 10⁻⁶ sterility assurance level (SAL):

| Cycle Type | Minimum Temperature | Minimum Exposure Time (Wrapped) | Minimum Dry Time |
|------------|---------------------|---------------------------------|-------------------------|
| Gravity | 132°C / 270°F | 15 minutes | 30 minutes ¹ |
| Pre-vacuum | 132°C / 270°F | 4 minutes | |

¹ Drying times vary according to load size and should be increased for larger loads.

7. Directions for use:

- A. Place the tip of the Chondral Pick on the surface of the bone at the point at which the socket/defect is to be created.
- B. As necessary use a suitable mallet to insert the Chondral Pick into the bone to the desired depth. The circumferential laser line indicates a distance of 3.5 mm from the tip.
- C. Remove the Chondral Pick

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Parcus Medical, LLC
MDSS GmbH