Anika Chondral Pick
Púa Condral Anika
Alêne chondrale
Anika Chondro Pick
Selezionatore condrale Anika
Cureta Condral Anika
Χονδροειδής Σφύρα Anika
„Anika“ chondralinis kirstuvas
Anika bruskspiss
Anika kondralhacka
Croșetă condrală Anika
Dłuto chrząstkowe Anika
Хондрален пикел Anika

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30990 Rev 5
Eff. 08/24

Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

ANIKA

www.anika.com
(941)755-7965

2797

EC REP

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30175 Hanover, Germany

Anika Chondral Pick (EN)

1. Intended Use:

The Anika Chondral Picks are designed to create a socket/defect in bone.

2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. Do not use this product if it is bent or there is visible damage to the shaft or drive section.
- C. This device is not intended to be used as an implant.
- D. Placing off-axis loads on the chondral pick may cause damage to the device.
- E. Caution: Common risks consistent with surgery may include: bone/tissue damage, revision surgery, soft tissue irritation, infection, allergic reactions.

3. Material:

This device is manufactured from a Stainless-Steel Alloy and may include an aluminum alloy handle. The materials used in the manufacture of this device are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Storage:

Products must be stored in a dry place.

6. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

7. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

² This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination

After sterilization aseptically transfer device into the sterile field in accordance with the institution's policy.

8. Directions for use:

- A. Aseptically transfer device into sterile field and place the tip of the Chondral Pick on the surface of the bone at the point at which the socket/defect is to be created.
- B. As necessary use a suitable mallet to insert the Chondral Pick into the bone to the desired depth. The circumferential laser line indicates a distance of 3.5 mm from the tip. Remove the Chondral Pick

9. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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