Draw Tight Suture Anchors
Anclajes de sutura Draw Tight
Ancres de suture Draw Tight
Draw Tight Fadenanker
Ancore da sutura Draw Tight
Âncora de Sutura Apertada de Sucção
Αγκίστρα ράμματος Draw Tight
„Draw Tight“ siūlų inkarai
Sammentrekbare suturankre
Draw Tight-suturankare
Ancore de sutură Draw Tight
Kotwy szewne Draw Tight
Фиксатори за стягане на конци

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

30999 Rev 7
Eff. 09/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 ANIKA  2797
www.anika.com
(941) 755-7965

 EC REP
MDSS GmbH
Schiffgraben 41
30175 Hanover, Germany

Draw Tight Suture Anchors (EN)

1. Device Description

The Draw Tight Suture Anchors are made with UHMWPE suture and a PEEK-OPTIMA power tip and are to be used for both large and small joint soft tissue repairs. The suture and anchor are provided preloaded onto a driver. The Draw Tight Suture Anchors are sterile and for single use only.

2. Intended Use

The Draw Tight Suture Anchors are indicated for attachment of soft tissue to bone.

3. Indications

The Draw Tight Suture Anchors are intended for the following indications:

<u>Shoulder</u>	Rotator Cuff Repair, Acromioclavicular Separation Repair, Bankart Lesion Repair, Biceps Tenodesis, Capsular Shift or Capsulolabral Reconstruction, Deltoid Repair, SLAP Lesion Repair.
<u>Knee</u>	Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Extra Capsular Reconstruction, Iliotibial Band Tenodesis, Patellar Ligament and Tendon Avulsion Repair.
<u>Foot/Ankle</u>	Lateral Stabilization, Medial Stabilization, Midfoot Reconstruction, Achilles Tendon Repair, Hallux Valgus Reconstruction, Metatarsal Ligament Repair.
<u>Elbow</u>	Tennis Elbow Repair, Biceps Tendon Reattachment.
<u>Hand/Wrist</u>	Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction, TFCC.

4. Contraindications

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the instructions for use or using the device for an indication other than those listed.

5. Intended Patient Population

The target patient groups for this device include patients requiring tissue repair involving the attachment of soft tissue to bone.

6. Intended Users

The intended users are healthcare professionals.

7. Intended Use Environment

The device is intended to be used in operating rooms in hospitals, outpatient surgical settings, and ambulatory surgery sites.

8. Clinical Benefits

The Draw Tight Suture Anchors aid in the fixation of soft tissue to bone during tissue repair, reattachment, and reconstruction procedures.

9. Device Lifetime

The functional lifetime of this device is 12-weeks which considers the time necessary for soft tissue healing.

10. Adverse Effects

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Common risks consistent with procedures using this type of device may include bone/tissue damage, intervention/revision surgery, soft tissue irritation, loss of fixation/repair, infection, allergic reaction(s).

11. Warnings

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.

- C. Anchor size selection should be made with care taking into consideration the quality of the bone in which the anchor is to be placed.
- D. Care and consideration should be given to the required depth of the prepared site and the anatomy of each patient in order to avoid damage to surrounding structures.
- E. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- F. Polyethylene polyblend sutures, as used in this device, elicit a minimal inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Polyethylene polyblend suture is not absorbed, nor is any significant change in tensile strength retention known to occur in vivo.
- G. The patient should be advised of the use and limitations of this device.
- H. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- I. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- J. This device must never be re-sterilized.
- K. Appropriate instrumentation should be used to implant this device.
- L. Loading the Draw Tight Suture Anchor in an off-axis manner may result in breakage or bending.

12. Precautions

- A. Excess pressure to distal tip of driver may damage driver shaft.
- B. Improper fixation of anchor may leave implant loose leading to soft tissue irritation, bone/tissue damage, and/or post operative failure.

13. Packaging and Labeling

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

14. Material Specifications

The Draw Tight Anchors are comprised of polyetheretherketone (PEEK) and high-strength, braided, polyethylene polyblend sutures. The Draw Tight Suture Anchors are MR-Safe.

15. Sterilization

The Draw Tight Anchor and driver assembly is supplied sterile. The contents are sterilized by EO gas. These products must never be re-sterilized.

16. Storage

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

17. Instructions for Use

- A. The Draw Tight Suture Anchor is supplied in a number of configurations. The associated instrumentation is keyed to the diameter of the suture anchor. See the table below for drill guide and drill bit selection.

Description	Reference #	Used With
8-Point Drill Guide (optional)	10446	Draw Tight 1.8mm Suture Anchor
1.8mm Draw Tight Drill Bit	10590	Draw Tight 1.8mm Suture Anchor
3.2mm Drill Guide (optional)	10591	Draw Tight 3.2mm Suture Anchor
3.2mm Draw Tight Drill Bit with Stop	10614	Draw Tight 3.2mm Suture Anchor
3.2mm Awl	10365	Draw Tight 3.2mm Suture Anchor

- B. Aseptically transfer device into sterile field and select the site for anchor placement ensuring that there is both the quantity and quality of bone to implant the anchor.
- C. Use one of the drill guides (as applicable) and drill bit /awl combinations from the table above to create a hole, in which to insert the anchor.

Note: When using a drill guide, the Draw Tight Suture Anchor is to be placed through the guide. The depth indicators on the driver, drill, and awl are keyed to the guides to insure proper anchor insertion depth.

- D. Taking care to follow the axis of the prepared hole, insert the Draw Tight Suture Anchor so that the circumferential line closest to the tip of the driver is flush with the surrounding bone.
- E. Remove the sutures from the driver handle and remove the driver and guide (if used).
- F. Pull the two ends of the white deployment suture firmly. This step deploys the anchor in the hole.
- G. Verify fixation by pulling on all sutures.
- H. The suture(s) marked with color will slide through the implant and are intended for fixation of tissue. Verify sutures slide freely by pulling on one end.
- I. If an additional suture is desired, the white deployment suture is available. If the white deployment suture is not used for fixation, cut the tails of the white deployment suture above the insertion site.
- J. Pass the desired sutures through the targeted tissue, fix the tissue with the surgeon preferred knots, and cut the suture tails above the knots.

18. Device / Instrumentation Disposal

Devices that have been removed from the packaging and not used or was explanted from a patient must be properly disposed of in accordance with the institution's policy.

19. Adverse / Reportable Event

Country	Authority	Email	Website
All events	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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