Non-cannulated Headed Reamers
Taladradoras con cabeza no canulada
Alésoirs dirigés, non tubulés
Nicht kanülierte Kopfreibahlen
Alesatori con testa non cannulati
Escareadores com Cabeça Não Canulada

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto

Instructions for Use
Instrucciones de uso
Mode d'emploi
Gebrauchsanweisung
Istruzioni per l'uso
Instruções de uso

31301 Rev 3
Eff. 09/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 **ANIKA**  **2797**
www.anika.com
(941) 755-7985

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Non-cannulated Headed Reamers (EN)

1. Intended Use

Non-cannulated Headed Reamers are used to create sockets and tunnels in bone.

2. Warnings

- A. Sterile devices designed for single use only should not be re-sterilized.
- B. Reuse of a device that is intended for single use could result in inferior device performance or patient infection.
- C. Never use a bent, dull, or damaged drill bit/reamer.
- D. Never use a device that has discoloration or staining.
- E. Never use a device with faded laser mark(s).
- F. Never use a damaged drill guide, guide pin or guide wire.
- G. These devices are not intended for use as an implant or to be left in the body.
- H. Applying excess pressure or using the device in an off-axis manner may result in breakage/bending or damage to both drill/reamer and guide pin/guide wire.
- I. Limit duration of use. Excessive use without irrigation may lead to heat buildup at the implant site.

3. Precautions

- A. Federal law restricts these devices to sale by or on the order of a physician.
- B. Stop using immediately if excessive force is required to advance the drill bit/reamer.
- C. Never use a device with any suspected defect or damage not covered in this IFU. If such damage is detected, notify the manufacturer.
- D. Common risks consistent with procedures using this type of device may include bone/tissue damage, intervention/ revision surgery, soft tissue irritation, loss of fixation/repair, infection, allergic reaction(s).

4. Material

This product is made from a stainless steel alloy. The material used in the manufacture of this device is radio-opaque, and can, therefore, be detected with conventional X-ray or Fluoroscopy.

5. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

6. Packaging and Labeling

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Contact the manufacturer to report any package/label damage or alterations.

7. Sterilization

These devices are sterilized by EO gas and supplied in sterile packaging. These devices must never be re-sterilized.

8. Storage

Devices must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Directions for Use

- A. Select a Non-cannulated Reamer of appropriate diameter.
- B. Open the product using aseptic technique.
- C. Place the non-fluted end of the reamer into an appropriately designed drill.
- D. Place the tip of the Non-Cannulated Reamer at the intended implant site.
- E. Prudently advance the Non-cannulated Reamer into the bone to the desired depth, using the drill motor set to rotate in a clockwise direction.

- F. Withdraw the Non-cannulated Reamer from the bone, remove it from the chuck of the drill and dispose of it in an appropriate manner.

10. Device Disposal

Sterile devices that have been removed from the packaging and reusable devices with visible damage(s) must be properly disposed of in accordance with the institution's policy.

11. Adverse / Reportable Events

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All events	Anika Therapeutics	Globalcomplaints@anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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