GFS BTB Product Family

Important Product Information

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GFS BTB Product Family (English)

1. Indications:

The GFS BTB devices are indicated for use in the fixation of ligaments and tendons in patients requiring ligament or tendon repair.

2. Contraindications:

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected, should be identified and precautions observed.
- D. Insufficient quality or quantity of bone.
- E. Patient's inability or unwillingness to follow surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the directions for use or using the device for an indication other than those listed.

3. Adverse Effects:

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Size selection of the implant should be made with care taking into consideration the quality of the bone on which the implant will rest and the desired length of graft that will reside in the tunnel.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. An adequate postoperative management plan should be implemented after implant removal.
- E. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations, are necessary components in achieving a good surgical result.
- F. This device, including all contents of the sterile package, must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- G. This device, including all contents of the sterile package, must never be re-sterilized.
- H. Appropriate instrumentation should be used to implant this device.
- I. The GFS BTB device must be used in conjunction with another GFS device, such as the GFS Mini or GFS Ultimate. While using the GFS BTB, it is important to ensure that all other GFS devices are being used as intended per their respective Instructions for Use (IFU).
- J. This device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment. Use of MR technology in the presence of devices of this nature may cause magnetically induced displacement forces and torques, radio frequency heating and image artifacts. Standard MRI screening guidelines for post-operative patients should be followed.

5. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications:

The GFS BTB Product Family is supplied with high-strength, braided, polyethylene polyblend sutures. The anchor material is Ti-6Al-4V ELI (ASTM F136).

7. Sterilization:



The GFS BTB Product Family is supplied sterile. The contents are sterilized by EO gas. These products must never be re-sterilized.

8. Storage:



Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Instructions for Use (Bone Tendon Bone Graft):

- A. Prepare bone tunnel using the user preferred method.
- B. Prepare the graft by drilling a hole in the bone block to be fixed into the femoral socket using a Parcus 3.0mm GFS BTB drill bit. Locate and center the hole 10mm from the proximal end of the bone block (opposite the graft), oriented in such a way that the hole passes through the cortical layer of the bone. It may be necessary to round the corners of the bone block to assist with orientation of the GFS BTB construct to the distal tip of the graft.
- C. Pass the GFS BTB through the bone block such that the button lies on the cortical face when the BTB suture is pulled through. Pass the BTB suture through the loop of the selected GFS Mini or GFS Ultimate device using the attached passing needle with suture.


- D. The passing needle is then used to pass the titanium button of the GFS BTB through the suture eye at the distal end of the GFS BTB suture construct. Once the GFS BTB is hooked back to itself, the passing needle may be cut off and discarded.

- E. Press the BTB button into the 3.0mm hole starting from the face which is opposite the cortical layer. From the cortical face of the bone block, pull the BTB suture using a non-abrasive instrument suitable for grasping suture until the button is fully inserted into the bone block. The

distal eyelet of the BTB suture should be pulling the button into the hole as this will prevent the suture from forming slack inside the hole.

- F. Once the GFS BTB is secured to the GFS Mini or GFS Ultimate, the construct is passed up through the previously created tibial and femoral bone tunnels following the recommendations of the applicable GFS device and surgeon technique preference.

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