## **Parcus Actiflip Product Family**

### **Important Product Information**

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## **Parcus Actiflip Product Family (English)**

### **1. Indications:**

The Parcus Actiflip is used for fixation of bone to bone or soft tissue to bone, and is intended as fixation post, a distribution bridge, or for distributing suture tension over areas of ligament or tendon repair in the knee, shoulder, and elbow and may include the following indications: anterior cruciate ligament, posterior cruciate ligament, pectoralis repair (minor/major), biceps tendon repair and reattachment (distal/proximal), acromioclavicular repair, and ulnar collateral ligament reconstruction.

The Parcus Actiflip is also intended for fixation of suture (soft tissue) to bone in the foot and ankle for midfoot reconstruction and hindfoot reconstruction with the following procedures: FHL Tendon Transfer, FDL Tendon Transfer, Posterior Tibialis Tendon Transfer, Anterior Tibialis Tendon Transfer.

### **2. Contraindications:**

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected, should be identified and precautions observed.
- D. Insufficient quality or quantity of bone.
- E. Patient's inability or unwillingness to follow surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the instructions for use or using the device for an indication other than those listed.
- G. The Actiflip Naked devices are intended to be used in conjunction with various types of suture. Surgeon selection of suture used with this device may impact the strength of the repair and could compromise fixation.

### **3. Adverse Effects:**

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

### **4. Warnings:**

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Size selection of the implant and suture should be made with care taking into consideration the quality of the bone on which the implant will rest and the anticipated forces that the fixation must resist.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. An adequate postoperative management plan should be implemented after implant removal.

- E. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- F. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- G. This device must never be re-sterilized.
- H. Parcus recommended instrumentation should be used to implant this device.
- I. This device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment. Use of MR technology in the presence of devices of this nature may cause magnetically induced displacement forces and torques, radio frequency heating and image artifacts. Standard MRI screening guidelines for post-operative patients should be followed.

**5. Packaging and Labeling:**

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

**6. Material Specifications:**

The Parcus Actiflip implant is comprised of Ti-6Al-4V ELI (ASTM F136). If applicable, supplied suture is high-strength, braided, polyethylene polyblend suture.

**7. Sterilization:**

This Parcus Actiflip Product Family is supplied sterile.

**8. Storage:**

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

**9. Instructions for Use:**

- A. A drill is required to create an appropriately sized hole in bone. Please refer to product packaging for recommended drill diameter.
- B. Based on the implant chosen and size of graft used, a stepped tunnel may be required.
- C. Using an appropriately sized Drill Tip Guide Pin or Drill Bit, create a bone tunnel to the necessary diameter and depth to accommodate the graft and fixation technique. A stepped tunnel may be required to insert the soft tissue graft into the bone socket. Surgeon selection, of the drill diameter is selected to closely match the graft diameter. A Parcus Graft Sizer or Parcus Graft Sizing Block can be used to assess graft diameter. Care must be exercised to preserve the distal cortex upon which the implant will rest.

**Note: If the lateral cortex of the bone is breached, use of a larger Actiflip**

**device may be required.**

- D. Attach the Actiflip implant to the graft. The free ends of the graft are whipstitched in the surgeon preferred manner, leaving adequate length of the suture tails that will be placed through the two center holes of the Actiflip implant. Aseptically open the appropriately sized Parcus Actiflip device and pass suture tails through center holes of the implant.
- E. Push the implant through the tunnel using the supplied insertion instrument. Apply light tension to the free sutures and press the button on the handle to release the implant from the insertion instrument. Pull tension on the sutures to ensure the implant is properly fixated in position.
- F. Complete the fixation by passing the sutures back through the graft and tying knots or by using adjunctive fixation, e.g. placing the sutures around a bicortical screw and washer.

**Suggested Technique for Subpectoral Biceps Tenodesis**

- A. Use the surgeon preferred surgical technique to expose the long head biceps tendon (LHBT) and identify the desired tenodesis site on the humerus.
- B. Disect and deliver the LHBT out of the incision.
- C. Using an appropriately sized drill bit, create a unicortical hole in the humeral cortex. A drill guide or tissue protector can be used to isolate the drill bit from the adjacent soft tissues.
- D. Insert the tip of the Actiflip into the drilled hole until it reaches the proximal aspect of the opposing cortical wall.
- E. Release button from Actiflip handle by depressing button on the handle.
- F. Leave the inserter in the hole and toggle/pull the attached sutures to aid in flipping the button in the medullary canal.
- G. Remove the inserter.
- H. Pull the LHBT through the double loops of suture attached to the Actiflip button and visualize desired length/tension point.
- I. Pass each suture end with needle through biceps tendon from deep to superficial once, with one proximal to the loops and the other distal to the loops.
- J. Pull on each free end of the suture together "cinching down" the loops/cradle until the LHBT is approximated firmly to the anterior humeral cortex.
- K. Tie standard knots.
- L. Cut and remove excess LHBT tissue.
- M. Irrigate wound and perform standard closure.

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