Parcus ATLAS

Important Product Information

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Parcus ATLAS (English)

1. Indications:

The Parcus ATLAS is indicated for CMC Arthroplasty.

2. Contraindications:

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected, should be identified and precautions observed.
- D. Insufficient quality or quantity of bone.
- E. Patient's inability or unwillingness to follow surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the instructions for use or using the device for an indication other than those listed.

3. Adverse Effects:

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings:

- A. Caution: Federal law restricts this device to sale by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. An adequate postoperative management plan should be implemented after implant removal.
- D. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- E. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- F. This device must never be re-sterilized.
- G. Parcus recommended instrumentation should be used to implant this device.
- H. This device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment. Use of MR technology in the presence of devices of this nature may cause magnetically induced displacement forces and torques, radio frequency heating and image artifacts. Standard MRI screening guidelines for post-operative patients should be followed.

5. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications:

The Parcus ATLAS implants are comprised of polyetheretherketone (PEEK), high-strength, braided, polyethylene polyblend suture and Ti-6Al-4V ELI (ASTM F136).

7. Sterilization:

Parcus ATLAS products are supplied sterile. The contents are sterilized by EO gas. These products must never be re-sterilized.

8. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Instructions for Use:

- A. A standard volar approach is used to expose the carpometacarpal (CMC) joint. Care is taken to identify and protect the branches of the superficial radial nerve.
- B. A longitudinal capsulotomy is made in the capsule, and radial flaps are elevated to expose the trapeziometacarpal joint.
- C. The trapezium is dissected and removed using the surgeon preferred technique.
- D. Aseptically open the ATLAS pouch and remove the Draw Tight 1.8mm Suture Anchor. Included within this package is the 1.8mm Draw Tight Drill Bit, GFS Micro Naked, and Suture Retriever. The base of the index metacarpal is identified and the 1.8mm Draw Tight Drill Bit is used to create an entrance into its medullary canal. The drill bit is advanced to the circumferential laser line above the distal end of the drill bit.
- E. The tip of the Draw Tight 1.8mm anchor is inserted into the previously drilled hole and advanced to the laser line on the driver shaft. Remove the sutures from the driver handle and remove the driver from the anchor.
- F. The white sutures attached to the anchor are identified. One of the white sutures is pulled until there is a solid resistance to the movement of the suture. The second white suture is pulled in the same fashion.
- G. Upon completing Step F, pull all of the sutures together to test fixation. Cut the white suture at the base of the second metacarpal.
- H. With care protecting soft-tissue structures, use the previously used 1.8mm Drill Bit to create a tunnel across the axis of the proximal base of the first metacarpal, drilling from radial to ulnar.
- I. Using the provided Suture Retriever, locate the Draw Tight sutures from the second metacarpal and pull them through the tunnel from ulnar to radial.
- J. Pass one suture through each of the center holes of the button. Tension the construct as desired, and then tie a secure knot across the button and cut suture tails.
- K. Close and dress the wound in the surgeon preferred manner.

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