Parcus Synd-EZ

Important Product Information

Wichtige Produktinformationen
Información Importante sobre el Producto
Informations Importantes sur le Produit
Informazioni Importanti sul Prodotto
Informações Importantes sobre o Produto
Informações Importantes sobre o Produto
Důležité informace o výrobku

Instructions for Use

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Parcus Synd-EZ (English)

1. Indications:

The Parcus Synd-EZ is intended to be used as an adjunct in fracture repair involving metaphyseal and periarticular small bone fragments where screws are not indicated, and as adjuncts in external and intramedullary fixation systems involving plates and rods, with fracture braces and casting. The Parcus Synd-EZ is intended to provide fixation during the healing process following syndesmotic trauma, such as fixation of the syndesmosis (syndesmosis disruptions) in connection with Weber B and Weber C ankle fractures.

2. Contraindications:

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected, should be identified and precautions observed.
- D. Insufficient quality or quantity of bone.
- E. Patient's inability or unwillingness to follow surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the instructions for use or using the device for an indication other than those listed.

3. Adverse Effects:

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings:

- A. Caution: Federal law restricts this device to sale by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and the potential of compromised results.
- C. The Synd-EZ device is intended to be used with or without a fracture plate. Special consideration to the necessity of a fracture plate and the quality of bone on which the implant will rest are imperative to the performance of this device.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. An adequate postoperative management plan should be implemented after implant removal.
- E. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations, are necessary components in achieving a good surgical result.
- F. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device and bone.
- G. This device, including all contents of the sterile package, must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices also causes cross-contamination leading to patient infection.
- H. Appropriate instrumentation should be used to implant this device.
- I. This device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment. Use of MR technology in the presence of devices of this nature may cause magnetically induced displacement forces and torques, radio frequency

heating and image artifacts. Standard MRI screening guidelines for post-operative patients should be followed.

- J. Special consideration should be made to the material compatibility of fracture plates and/or screws that may be used in conjunction with the Synd-EZ device. Parcus Medical recommends that implants in direct contact with each other should be comprised of identical metallurgical composition.

5. Packaging and Labeling:

- A. Do not use this product, if the packaging or labeling has been damaged or altered in any way, or if there is any evidence of exposure to moisture or extreme temperature.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications:

The Synd-EZ is supplied with high-strength, braided, polyethylene polyblend sutures.

- A. The Synd-EZ Ti oblong anchor, top hat button, and washer are made of Ti-6Al-4V ELI (ASTM F136).
- B. The Synd-EZ SS top hat button material is 316LS/316LVM (ASTM F138) and the oblong anchor is made of Ti-6Al-4V ELI (ASTM F136).

7. Sterilization:

The Synd-EZ is supplied sterile. The contents are sterilized by EO gas. These products must never be re-sterilized.

8. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Instructions for Use (Syndesmosis Repair):

- A. Make an incision on the lateral side of the ankle. A fracture plate may be required to stabilize and realign the fibula.
- B. Aseptically open the Synd-EZ package. If a fracture plate is necessary, select a hole in the plate per surgeon preference and drill a 1.5mm K-Wire through the center of the selected hole. If a fracture plate is not required, drill a 1.5mm K-Wire in the surgeons preferred location. The K-Wire is advanced from lateral to medial until it protrudes through the skin on the medial side.
- C. Leaving the K-Wire in position, advance the 3.5mm cannulated drill bit from lateral to medial through the medial cortex of the tibia.
- D. Remove the K-Wire and drill bit.
- E. If using the Synd-EZ Ti without a fracture plate, the provided Titanium washer may be installed at this time to increase the surface area of the top hat button. It is noted that the Synd-EZ SS is not provided with a washer. To install the washer, thread the solid white Passing suture-tail, white/black Toggle suture tail, and oblong anchor through the countersunk surface of the washer. (**NOTE:** It is imperative that the washer is installed such that the top hat button sits flush inside of the countersink).
- F. Thread the solid white Passing suture-tail and the white/black Toggle suture-tail through the eyelet in the 2.4mm passing pin.
- G. Insert the tip of the passing pin into the drilled hole. Carefully advance the passing pin from lateral to medial until the oblong button of the Synd-EZ construct is resting on the medial side of the tibia. Remove suture tails from the eyelet and discard the passing pin.
- H. Use the white/black suture loop (Toggle suture) and solid white suture loop (Passing suture) to position the oblong button so that it is perpendicular to the 3.5mm tunnel and resting on the medial cortical surface of the tibia.

- I. On the top hat button, pull the solid blue Counter-traction suture loop to tension the construct and verify that the oblong button has turned perpendicular to span the hole.
- J. Place the white/blue Adjustment suture loop into the central groove in the black suture pull handle.
- K. Pull the black suture pull handle with consistent force to slowly repel the top hat button into the hole in the fracture plate (or bone). (**NOTE:** In intervals of approximately 5-10mm, lightly pump the solid blue Counter-traction suture loop to maintain the alignment of the suture construct and top hat button. This will encourage a steady, controlled sliding motion of the top hat. Do not simultaneously pull the Adjustment and Counter-traction suture loops as this will increase the amount of force needed to repel the top hat button.)
- L. When the Adjustment suture is tight, the top hat button should be seated flush or slightly proud of the fracture plate or cortical surface of the fibula. The suture construct should be tight. Confirm fixation by tensioning the solid blue Counter-traction suture loop while inspecting the top hat for excess movement. Apply a final, firm pull on the black suture pull handle.
- M. Complete the syndesmotic repair by cutting and removing the Toggle/Passing/Counter-traction/Adjustment suture loops.
- N. If desired, an additional Synd-EZ may be used to reinforce fixation or alignment of the tibiofibular joint.

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