## Parcus Suture Pull Handle

### Important Product Information

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## Parcus Suture Pull Handle (English)

### 1. Indications:

The Parcus Suture Pull Handle is designed as an ancillary instrument to assist in the tensioning and handling of Parcus products with an adjustable suture.

### 2. Warnings:

- A. Caution: Federal law restricts this device to sale by or on the order of a physician.
- B. Parcus Suture Pull Handles are supplied unsterile as reusable instruments.
- C. Never use a deformed or damaged suture pull handle.
- D. This product is not intended for use as an implant.
- E. Be sure to remove any debris from the cannulation of the handle during the cleaning process.

### 3. Material:

This instrument is manufactured from Delrin/Acetal, is radiolucent, and may be visualized through the use of ultrasound or MRI. No part of this device is intended to be implanted in the body.

### 4. Inspection:

- A. Inspect the instrument for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

### 5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub instrument with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.
- D. The instrument should be placed in an ultrasonic cleaning unit and processed according to the unit's directions.
- E. The instrument should be rinsed thoroughly with water following the ultrasonic process.

### 6. Sterilization:

This instrument is not provided sterile. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a  $10^{-6}$  sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 minutes	30 minutes <sup>1</sup>
Pre-vacuum	132°C / 270°F	4 minutes	

<sup>1</sup> Drying times vary according to load size and should be increased for larger loads.

### 7. Instructions for Use:

- A. Create a suture loop using the surgeons preferred method in the adjustable suture of a suture anchor or suture suspension device. Note: Some anchors and suture suspension devices are supplied with a suture loop and may not require an additional suture loop.
- B. Locate the central circumferential groove in the suture pull handle.
- C. Take the loop in one hand and the suture pull handle in the other and slide the suture loop over-top of the suture pull handle.
- D. Place the loop such that it is resting in the central circumferential groove in the suture pull handle.
- E. Pull tension on the handle as necessary and follow the IFU that is provided with the applicable anchor or suture suspension device.

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