Parcus PTS Suture Passing Needle

Important Product Information

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Parcus PTS Suture Passing Needle (English)

1. Indications:

The Parcus PTS Suture Passing Needle is designed to be used with the Parcus PTS Suture Passing Hand Piece to affect the placement of #2 suture through soft tissue.

2. Warnings:

- A. Caution: Federal law restricts this device to sale by or on the order of a physician.
- B. This product is designed for single use only.
- C. Never use a PTS Suture Passing Needle that has surface defects or if the plastic hub is damaged in any way. PTS Suture Passing Needles that are reused can fatigue and break.
- D. This product is not intended for use as an implant.
- E. The Parcus PTS Suture Passing Needle is used in conjunction with the Parcus PTS Suture Passing Hand Piece, which is supplied separately.
- F. Ensure that the PTS Suture Passing Needle is retracted into the PTS Suture Passing Hand Piece before passing it into the body. Use of an appropriately sized cannula (not supplied) is recommended to facilitate insertion of the device. A damaged needle can result in product failure.

3. Material:

This device is manufactured from a Nitinol Alloy with a molded plastic hub. The materials used in the manufacture of this device that are intended to be placed inside the body, are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature, or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

5. Sterilization:

The Parcus PTS Suture Passing Needle is supplied sterile.

6. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

7. Instructions for use:

- A. Aseptically open the packaging containing the Parcus PTS Suture Passing Needle (Needle).
- B. Place the Needle into a sterile Parcus PTS Suture Passing Hand Piece (Hand Piece) by aligning the shaft of the Needle with the luer lock located on the Hand Piece.
- C. Insert the tip of the Needle into the cavity of the luer lock and slide it forward until the plastic hub of the Needle can be engaged in and secured in the slots of the Hand Piece.
- D. Actuate the Needle by grasping and squeezing the trigger of the Hand Piece. The Needle should pass out and up through the capture window of the upper jaw of the Hand Piece and retract back into the Hand Piece when the trigger is released. **Do not use the device if the Needle fails to function in this manner.** If the Needle does not function properly, replace the Needle and repeat the test. If it fails with a new Needle, check the Hand Piece for damage. Do not use a damaged Hand Piece.
- E. Load suture or suture tape into the slot located on the right side of the lower jaw of the Hand Piece. The suture tail should be oriented on the bottom of the lower jaw.
- F. Be sure to squeeze the trigger on the Hand Piece to close the jaw, but not enough to expose the Needle, and insert the Hand Piece into the joint through an appropriately sized cannula.

- G. Once in the joint, open the jaw by releasing the trigger and position the targeted tissue between the upper and lower jaws of the Hand Piece.
- H. Squeeze the trigger of the Hand Piece with enough force to grasp the tissue and deploy the Needle.
 - i. If using the self-capture configuration of the device, the Needle will push the suture into the upper jaw of the Hand Piece.
 - ii. When using non-capture device, the instrument may be removed, leaving the suture proud in the targeted tissue.
- I. To release the captured suture, pull the body of the suture towards the trigger of the Hand Piece.

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