Parcus Tap without Handle

Important Product Information

Produktinformationsblatt

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Parcus Tap without Handle (English)

1. Indications:

The Parcus Tap without Handle is designed to prepare bone for the implantation of the Parcus Suture Anchors. The shafts are designed to fit into handles/means of insertion that utilize an AO chuck.

2. Warnings:

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. This product should not be used to implant any device other than the Parcus Suture Anchors.
- C. This product is not intended for use as an implant.
- D. Do not use this product if it is bent or there is visible damage to the tip or thread form.

3. Material:

This device is manufactured from a Stainless Steel Alloy shaft. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

6. Sterilization:

- A. This device is not provided sterile. It should be cleaned and sterilized using validated protocols before each use.
- B. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a 10^{-6} sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

7. Instructions for use:

- A. Connect the shaft of the Parcus Tap to the intended handle/means of insertion, ensuring that it is securely seated.
- B. Identify the site for anchor insertion and prepare using either a drill or an awl per the technique recommended for the Parcus Suture Anchor to be used. Press the Parcus Tap into the bone until the threads can be engaged.
- C. Rotate the Tap clockwise into the bone until the top of the laser line is flush with the bone surface.
- D. Remove the Tap by rotating in a counter-clockwise direction until its threaded portion is clear of the bone.
- E. Insert a Parcus Suture Anchor into the prepared hole per the technique identified in the Instructions for Use provided with the Suture Anchor.

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