Anika Twist Biocomposite Tap
Terraja de biocompuesto Anika Twist
Taraud biocomposite Anika Twist
Anika Twist Biocomposite Gewindebohrer
Maschio per biocomposito Anika Twist
Torneira de Biocompósito de Torção Anika
Βιοσύνθετο Βύσμα Twist Anika
„Anika Twist“ biokompozito sriegiklis
Anika spiralgjengetapp av biokompositt
Anika vridbar biokomposittapp
Burghiu Anika pentru Twist Biocomposite
Wkręcany gwintownik biokompozytowy Anika
Биокomпозитен метчик Anika Twist

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Instructions for Use
Instrucciones de uso
Mode d'emploi
Gebrauchsanweisung
Istruzioni per l'uso
Instruções de uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Användarinstruktioner
Instrucțiuni de utilizare
Instrukcja użytkowania
Инструкции за употреба

31798 Rev 5
Eff. 08/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 **ANIKA**
www.anika.com
(941)755-7965

 2797

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Anika Twist Biocomposite Tap (EN)

1. Intended Use:

The Anika Twist Biocomposite Tap is designed to prepare bone for the implantation of the Anika Twist Biocomposite Suture Anchors.

2. Warnings:

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. This product should not be used to implant any device other than the Anika Suture Anchors.
- C. This product is not intended for use as an implant.
- D. Do not use the Tap if it is bent or there is visible damage to the tip or threads.

3. Material:

This device is manufactured from a Stainless-Steel Alloy shaft and an Aluminum handle. The materials are radio-opaque and can be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Completely submerge used instruments in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the instrument until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush)
- B. Remove instrument from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge instrument in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the instrument or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the instrument is free of blood or soil before moving to step F.
- F. Remove excess moisture from the instrument with a clean, absorbent, and non-shedding wipe.

6. Sterilization:

This device is not provided sterile. Devices shall be double wrapped in wrapping material or pouches, double wrapped while in a container or tray, or placed within a reusable sterilization container. Combining two or more systems may hinder the efficacy of the sterilization to an unacceptable level. Reusable devices shall be sterilized with a load size which does not exceed 10,630 grams. The following table provides recommended minimum sterilization parameters that have been validated by Anika Therapeutics to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹
	134°C / 273°F	3 minutes	
	132°C / 270°F	18 minutes ²	50 minutes ¹

¹Drying times vary according to load size and should be increased for larger loads.

²This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern regarding TSE/CDJ contamination. This cycle is not to be used for the inactivation of TSE/CJD contamination

7. Instructions for use:

- A. Identify the site for anchor insertion and press the Anika Twist Biocomposite Tap into the bone until the threads can be engaged. Use a suitable mallet if necessary.

- B. **Note:** If desired, the quality of the patient's bone may be evaluated with the use of a Anika Awl prior to tapping. The 4.5mm Twist Biocomposite Soft Bone Tap requires the use of an awl or drill bit prior to tapping. Reference the IFU provided with the suture anchor for the appropriate sizing.
- C. Rotate the Tap clockwise into the bone until the top of the laser line is flush with the bone surface.
- D. Remove the Tap by rotating the handle in a counter-clockwise direction until its threaded portion is clear of the bone.
- E. Insert a Anika Suture Anchor into the prepared hole.

8. Adverse/ Reportable Event:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All	Anika Therapeutics	ParcusCS@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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