ProPass Suture Passer
Pasador de suturas ProPass
Passe-suture ProPass
ProPass Fadentransporter
Passasutura ProPass
Passador de Sutura ProPass
Συσκευή Διέλευσης Ραμμάτων ProPass
„ProPass“ siūlų tempiklis
ProPass suturpasser
ProPass suturförare
Dispozitiv de trecere a suturii ProPass
Przeciągacz do nici ProPass
Уред за прокараване на шевове ProPass

Important Product Information
Información importante sobre el producto
Informations importantes relatives au produit
Wichtige Produktinformationen
Informazioni importanti sul prodotto
Informações Importantes sobre o Produto
Σημαντικές πληροφορίες για το προϊόν
Svarbi informacija apie gaminį
Viktig produktinformasjon
Viktig produktinformation
Informații referitoare la produs
Ważne informacje o produkcie
Важна информация за продукта

Directions for Use
Instrucciones de uso
Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instruções para o uso
Οδηγίες χρήσης
Naudojimo instrukcija
Bruksanvisning
Anvisningar för användning
Instrucțiuni de utilizare
Wskazówki dotyczące użytkowania
Инструкции за употреба

31960 Rev 3
Eff. 09/24

 Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243 USA

 ANIKA  2797
www.anika.com
(941) 755-7965

 EC REP
MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

ProPass Suture Passer (EN)

1. Device Description:

The ProPass Suture Passer includes a reusable suture passer (handpiece) and a disposable suture passer needle which are sold separately.

The ProPass Suture Passer is a handheld stainless-steel device consisting of a 2-piece spring loaded handle, a cannulated shaft, tissue grasping jaws, and a spring-loaded trigger that activates or releases the grasper's upper jaw. The Suture Passer is available with or without a spring-loaded trap door on the upper jaw to capture suture after deployment. The ProPass Suture Passer is used in conjunction with the ProPass Suture Passer Needle to pass and deploy #2 braid suture and 1.6mm suture tapes.

The ProPass Suture Passer Needle is a disposable device comprised of a sharp and flexible notched needle ribbon and a plastic hub that fits into the spring-loaded handle of the ProPass Suture Passer. The suture passer needle is a single-use device that is provided sterile.

2. Intended Purpose/Indication:

The ProPass Suture Passer is intended to affect the passage and deployment of suture through soft tissue.

3. Contraindications:

None.

4. Device Performance/Lifetime:

The ProPass Suture Passer handpiece is a reusable device and may be used if no damage is visible/evident. The ProPass Suture Passer Needle is a single use, disposable device.

5. Warnings:

ProPass Suture Passer Needle

- A. This device is designed for single use only.
- B. This device must never be re-used. Suture passer needles that are reused can fatigue and break. Reuse or reprocessing of single use devices can cause cross-contamination leading to patient infection.
- C. Never use a suture passer needle that has surface defects or if the plastic hub is damaged in any way. A damaged needle can result in device failure.
- D. This device is not intended for use as an implant or to be left in the body.
- E. Ensure the suture passer needle is retracted into the suture passer before approaching the surgical site.
- F. This device must never be re-sterilized.

ProPass Suture Passer

- G. Always use the suture passer under direct or arthroscopic visualization.
- H. Do not use if the device is damaged.
- I. Do not use the suture passer on bone or other similarly hard tissue, as this may cause wear to the device or patient injury.
- J. Avoid activating the needle multiple times without tissue inside the grasping jaw, to minimize potential for wear or damage to the needle.
- K. The upper jaw should be activated to firmly grasp and fully reduce the soft tissue against the lower jaw before advancing the needle. Failure to firmly grasp and reduce the soft tissue upon needle activation could result in needle breakage.
- L. Do not apply excessive axial or bending forces during use. Excessive force may limit the functionality or cause the handpiece and/or needle to bend, deform, or break.

6. Cautions:

- A. Federal law restricts this device to be sold by or on the order of a physician.
- B. Do not use on patients with known allergies and other reactions to device materials.
- C. Common risks consistent with this type of device may include bone/tissue damage, revision surgery, soft tissue irritation, loss of fixation/repair, infection, allergic reaction(s).

ProPass Suture Passer Needle

- D. Never use a device with any suspected defect or damage. If such damage is detected, notify the manufacturer.
- E. The device has sharp edges that could tear and puncture surgical gloves.
- F. The device has sharp edges that could cause unintended damage. Use caution around the articular surface and surrounding tissues.
- G. The needle should be retracted immediately when removing the device. Failure to do so could result in device damage and/or patient injury.

ProPass Suture Passer

- H. Inspect the device and packaging prior to use for any visible damage and proper mechanical function.
- I. Ensure proper function by opening and closing the jaws or deploying and retracting the needle one time. If damage is detected, do not use the device and consult the manufacturer for guidance.
- J. A damaged or broken device may result in extended or additional surgeries or residual foreign bodies.
- K. Ensure tissue grasping jaws have proper clearance before the needle is activated. Do not activate the jaws inside an arthroscopic cannula.
- L. Improper activation of the suture passer may fail to pass/deploy suture into soft tissue.
- M. If using a cannula, check for free passage of the device. It is recommended to always use a cannula (7mm diameter or larger) in conjunction with the suture passer.

7. Material Specifications:

The ProPass Suture Passer is manufactured from a Stainless-steel alloy. The ProPass Suture Passer Needle is manufactured from Nitinol with a machined plastic hub. The materials used in the manufacture of these devices are radio-opaque and can be detected with conventional X-Ray or fluoroscopy.

Hazardous substances (If applicable): None

8. Packaging and Labeling:

- A. Do not use this device if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature, or has been altered in any way.
- B. Please contact Anika Customer Service to report any package/label damage or alterations.

9. Sterilization:

ProPass Suture Passer

- A. The ProPass Suture Passer is provided non-sterile and should be cleaned and sterilized using validated protocols before each use.
- B. The following table provides the recommended minimum sterilization parameters to provide a 10^{-6} Sterility Assurance Level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time*
Pre-Vacuum	132°C	4 minutes	30 minutes
	134°C	3 minutes	
(WHO) Steam Cycle^	134°C	18 minutes	50 minutes

*Drying times vary according to load size and should be increased for larger loads.

^This cycle is not for use in the United States. These are disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments/devices where there is a concern regarding TSE/CJD contamination. This cycle is not to be used for the inactivation of TSE/CJD contaminants.

ProPass Suture Passer Needle

The ProPass Suture Passer Needle is supplied in sterile packaging. The contents are sterilized by ethylene oxide (EO) gas.

10. Storage:

Devices must be stored in the original unopened package in a dry place. The ProPass Suture Passer Needle must not be used beyond the expiration date indicated on the package.

11. General Instructions for use:

Insert Needle

1. Aseptically open the packaging containing the ProPass Suture Passer Needle. With the thicker portion of the needle hub facing up, grasp the shaft of the needle and introduce the distal tip into the proximal aperture of the suture passer cannulated shaft.
2. Advance the needle until the hub settles into the proximal slot on the superior portion of the handle and the hub's inferior recess engages the cross pin within the slot.

Load Suture

3. Keep jaws open for suture loading. Suture is loaded on the right-hand side directly into the front of the bottom jaw, with suture tail down or tail up. Pull the suture back proximally to seat it into the slot recess.
CAUTION: Use only standard arthroscopic surgical suture and ProPass Suture Passer Needle. If using the device in an open surgery, wet the distal end of the device with sterile water.

Grasping Tissue

4. With the suture loaded, the jaw of the suture passer closed, and the needle retracted in its resting position, pass the suture passer into the surgical site under arthroscopic or direct visualization.
CAUTION: Hold the suture slightly taut to ensure that it does not become entangled on itself or in the suture passer while it is passed into the surgical site.
5. Open the jaw of the suture passer, grasp and hold the targeted tissue for suture passing by pulling the FRONT trigger (jaws movement). Close the jaw over the tissue with a FIRM GRASP.

Fire Needle

6. Hold the front trigger in the closed position and squeeze the device handle, advancing the needle through the soft tissue to deliver suture through the tissue.
7. Release squeeze on the handpiece to allow the needle to deploy the suture on the adjacent side of the targeted tissue and to return the needle to its resting position inside the handpiece.
NOTE: If significant force is required to pass the needle, loosen the grip on the tissue and try again. If a false start occurs and the needle is partially deployed, the device may need to be removed and the suture reloaded.

Retrieve and Remove Suture

8. If using the non-auto capture configuration, retrieve the deployed suture with a suture grasper.
9. Release the grasped tissue and slowly remove the ProPass Suture Passer from the surgical site.
CAUTION: If the needle is not fully retracted, the passer will be difficult to remove from the tissue.
10. Repeat steps 4-7 for each suture passing.
11. If using the auto capture configuration, pull the captured suture proximally and inferiorly (in relation to the device).
12. If suture does not release from the trap door, slightly squeeze the handle to advance the needle enough to release the trap door's grasp while simultaneously pulling suture limbs proximally and inferiorly (in relation to the device).
13. When the preparation is finished, remove and dispose of the needle.

12. Inspection:

- A. Devices intended to be reused should be inspected prior to each use to assess damage, wear, and to ensure all components are functioning as intended. Signs of damage and wear include surface discoloration, surface corrosion, surface damage, general bending, cracking, and fractures.
- B. Some inspection points include:
 - a. Jaw actuation via the trigger;
 - b. The suture passer cannulation;
 - c. Actuation of the suture passer needle, with jaw closed via the handle.
- C. If the device does not pass inspection, do not use, and notify an Anika representative.

13. Reuse/Cleaning Instructions:

ProPass Suture Passer handpiece

- A. Completely submerge used devices in enzyme solution and allow to soak for 20 minutes. Use a soft bristled, nylon brush to gently scrub the device until all visible soil has been removed. Particular attention must be given to crevices, lumen, mated surfaces, connections and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft bristled brush (i.e. pipe cleaner brush).
- B. Remove device from enzymatic solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- C. Place prepared cleaning agents in a sonication unit. Completely submerge device in cleaning solution and sonicate for 10 minutes at 45-50 kHz.
- D. Rinse device in purified water for at least 3 minutes or until there is no sign of blood or soil on the device or in the rinse stream. Thoroughly and aggressively flush lumens, holes, and other difficult-to-reach areas.
- E. Repeat steps A through D a second time and ensure the device is free of blood or soil before moving to step F.
- F. Remove excess moisture from the device with a clean, absorbent, and non-shedding wipe.

ProPass Suture Passer Needle: Not applicable

14. Device Disposal:

Single use devices that have been removed from the packaging and reusable devices with visible damage must be properly disposed of in accordance with the institution's policy.

15. Adverse/Reportable Events:

Any serious incident that occurs in relation to the device should be reported to both Anika Therapeutics and the applicable authority using the following contact information.

Country	Authority	Email	Website
All events	Anika Therapeutics	Globalcomplaints@Anika.com	www.anika.com
Australia	Therapeutic Goods Administration (TGA)	IRIS@health.gov.au	www.tga.gov.au/reporting-adverse-events
European Union	EC Representative	info@mdssar.com	www.mdssar.com

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