
Instructions for Use

ORTHOVISC®-T Sodium Hyaluronate

FOR PERIOSTEOTENDINOUS INJECTION

This instructions for use is intended exclusively for distribution outside of the USA.

MANUFACTURER CONTACT INFORMATION

Anika Therapeutics, Inc.
32 Wiggins Avenue
Bedford, MA 01730
USA

AUTHORIZED EU REPRESENTATIVE CONTACT INFORMATION

Anika Therapeutics S.r.l.
Corso Stati Uniti, 4/U
35127 Padova (PD)
ITALY

In an effort to best meet the needs of our customers and minimize waste, these Instructions for Use are being provided in electronic format. This document is subject to change; the most current version of this IFU is available online. If unsure if using the latest revision, please reprint the IFU at www.anikaifu.com. If a paper copy is preferred it may be requested, free of charge, by contacting Anika Therapeutics, Inc. at www.anikaifu.com. The onus resides with the user to ensure that the most up-to-date IFU is used.

Products may not be available in all markets. Product availability is subject to the regulatory clearance in individual markets. Please contact your local representative if you have questions about product availability in your area.

ORTHOVISC-T PACKAGE INSERT**DESCRIPTION**

Orthovisc®-T is a sterile viscoelastic preparation supplied in a disposable glass syringe containing 2.0 ml of sodium hyaluronate dissolved in physiological saline. The hyaluronan is derived from bacterial fermentation and contains no animal material.

CHARACTERISTICS

Sodium hyaluronate is a high molecular weight polysaccharide composed of sodium glucuronate and N-acetylglucosamine. Sodium hyaluronate is ubiquitously distributed throughout the tissues of the body and is present in high concentrations in synovial fluid and the extracellular matrix. Sodium hyaluronate functions as a tissue lubricant. The sodium hyaluronate in Orthovisc®-T promotes tendon gliding and provides an environment to support tendon repair. Orthovisc®-T is biocompatible and non-inflammatory.

INDICATIONS

The Orthovisc®-T device is intended to relieve pain and restore function in tendons affected by chronic lateral epicondylitis. It is not intended to be used for the treatment of first or second degree ankle sprain or symptomatic rotator cuff tendinopathy.

CONTRAINDICATIONS

Orthovisc®-T should not be used in patients who have known hypersensitivity to any of its constituents.

DOSAGE AND ADMINISTRATION

Orthovisc®-T is a medical device that should be used by a physician trained in periosteotendinous injections only. Orthovisc®-T should be injected around the affected osteotendinous junction once a week for a total of two injections. Injections may be done under imaging guidance.

The contents of the Orthovisc®-T syringe are sterile as long as the Orthovisc®-T pouch remains unopened and undamaged. Remove the pre-filled syringe from the pouch, unscrew the Luer-lock cap, attach a suitable disposable needle (e.g. 22 – 25 gauge), and secure by turning slightly. Discard the syringe and needle after single use.

PRECAUTIONS

General precautions for periosteotendinous injections should be observed. Strict aseptic injection technique should be used during the application of Orthovisc®-T. Avoid injections into nerves or blood vessels. Do not use Orthovisc®-T if the package has been opened or damaged. The pre-filled syringe is intended for single use only. The contents of the syringe should be used immediately after opening. Discard any unused Orthovisc®-T. Do not resterilize.

ADVERSE REACTIONS

Sodium hyaluronate is a natural component of the tissues of the body and is non-inflammatory; therefore any inflammatory response is considered to be caused by the injection procedure itself. There are minimal risks associated with periosteotendinous injections, including injection site pain, swelling, erythema, bleeding, and infection.

HOW SUPPLIED

Orthovisc®-T is a sterile viscoelastic preparation of sodium hyaluronate supplied in a disposable glass syringe. Each mL of Orthovisc®-T contains 15 mg of sodium hyaluronate, 9 mg of sodium chloride and sterile water for injection (USP). Orthovisc®-T is sterile filtered and aseptically filled. The content of the syringe is processed in an aseptic manufacturing environment and are sterile to an SAL level of 10^{-3} if the syringe is intact. The Orthovisc®-T syringe is packaged in a pouch. Do not use Orthovisc®-T if the pouch has been opened or damaged. Orthovisc®-T should be stored at 2°C to 25°C, and should be allowed to reach room temperature approximately 20-45 minutes prior to use. Protect from freezing.

SHELF LIFE

The shelf life of Orthovisc®-T is 24 months when stored at a temperature of 2°-25°C.

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