en

Instructions for Use

MONOVISC[®]

Lightly Cross-linked High Molecular Weight Sodium Hyaluronate

This instructions for use is intended exclusively for distribution outside of the USA.

In an effort to best meet the needs of our customers and minimize waste, these Instructions for Use are being provided in electronic format. This document is subject to change; the most current version of this IFU is available online. If unsure if using the latest revision, please reprint the IFU at www.anikaifu.com. If a paper copy is preferred it may be requested, free of charge, by contacting Anika Therapeutics, Inc. at www.anikaifu.com. The onus resides with the user to ensure that the most up-to-date IFU is used.

Products may not be available in all markets. Product availability is subject to the regulatory clearance in individual markets. Please contact your local representative if you have questions about product availability in your area.

Read this instruction completely prior to use.

DESCRIPTION

MONOVISC® is a sterile, non-pyrogenic, sodium hyaluronate solution, lightly cross-linked with a proprietary chemical cross-linker. MONOVISC contains 22 mg/mL of lightly cross-linked sodium hyaluronate (NaHA) dissolved in phosphate buffered saline at physiological osmolality. MONOVISC is manufactured from ultra pure, high molecular weight sodium hyaluronate produced by bacterial fermentation. Hyaluronic acid is a natural complex polysaccharide of the glycosaminoglycan family.

CHARACTERISTICS

Sodium Hyaluronate is a high molecular weight polysaccharide composed of sodium glucuronate and N-acetylglucosamine. The sodium hyaluronate in MONOVISC is derived from bacterial fermentation. Hyaluronic acid is ubiquitously distributed throughout the tissue of the body and is present in high concentrations in such tissues as vitreous humor, synovial fluid, umbilical cord and dermis. Sodium hyaluronate functions as a tissue lubricant and is thought to play an important role in modulating the interactions between adjacent tissues. It can also act as a viscoelastic support maintaining a separation between tissues. Different sodium hyaluronate preparations may have different molecular weights but have the same chemical structure. The MONOVISC cross-linked sodium hyaluronate injection is biocompatible, non-inflammatory and non-pyrogenic. Sodium hyaluronate preparations have been shown to be well tolerated in osteoarthritic synovial joints such as the knee, hip, shoulder, and ankle. The actions of MONOVISC are lubrication and mechanical support.

APPLICATION

MONOVISC Cross-Linked Sodium Hyaluronate Injection is a single, intra-articular injection of cross-linked sodium hyaluronate designed to treat the symptoms of osteoarthritis.

INDICATIONS

MONOVISC is indicated as a viscoelastic supplement or a replacement for synovial fluid in human knee, hip, shoulder, and ankle joints. MONOVISC is well suited for treatment of the symptoms of mild to moderate osteoarthritis of the knee, hip, shoulder, and ankle joints for patients who have failed to respond adequately to conservative non-pharmacologic therapy and simple analgesics. In clinical studies, Monovisc has been proven to have a duration of effect of at least 6 months and a second injection 6 months after the first injection was shown to be safe.

DIRECTIONS FOR USE

The required amount of MONOVISC is injected through a sterile, disposable, hypodermic needle of suitable gauge into the selected joint space. The sterile needle should be attached to the MONOVISC syringe by a health care professional using a health care facility-approved aseptic technique. Common needle gauges for injections are 18-21 gauge. The final needle selection for any procedure is determined by the physician. The health care provider should ensure proper penetration into joint synovial space prior to injecting MONOVISC.

CONTRAINDICATIONS

MONOVISC is composed of cross-linked sodium hyaluronate and may contain trace amounts of gram positive bacterial proteins. The following pre-existing conditions may constitute relative or absolute contraindications to the use of MONOVISC:

- Known sensitivity to any of the materials contained in MONOVISC
- Pre-existing infections of the skin region of the intended injection site
- Known infection of the index joint
- Known systemic bleeding disorders
- Weak populations, such as children, pregnant or lactating women

PRECAUTIONS

- Those precautions normally considered during injection of substances into joints are recommended.
- Only medical professionals trained in accepted injection techniques for delivering agents to joint spaces should inject sodium hyaluronate for this application.

- The amount of MONOVISC necessary to be injected depends on patient anatomy and needs to be defined by the medical professional performing the procedure. An excess quantity of sodium hyaluronate should not be used and the patient should be monitored closely.
- The synovial space should not be overfilled.
- If pain increases during the injection procedure, the injection should be stopped, and the needle withdrawn.
- Transient pain or swelling may occur after the intra-articular (IA) injection.
- As with any invasive joint procedure, it is recommended that patients avoid strenuous or prolonged (i.e., more than one hour) weight-bearing activities such as running or tennis within 48 hours following the intra-articular injection.

SSCP

The summary of the safety and clinical performance (SSCP) is available in EUDAMED (<https://ec.europa.eu/tools/eudamed>) which is linked to the product Basic UDI- DI, 00817337MONOVGU.

ADVERSE REACTIONS

Adverse reactions to intra-articular injection of Monovisc based on clinical study data tend to be transient in nature and include the following:

•Joint Pain/Arthralgia •Injection Site Pain •Swelling •Stiffness •Rash •Effusion •Injection Site Reaction •Redness
•Pain/Pain in Extremity •Procedural Pain •Bruising •Synovitis •Subcutaneous Nodule •Cyst

HOW SUPPLIED


MONOVISC is a sterile viscoelastic preparation supplied in a disposable glass syringe delivering 4.0 mL. Each mL of MONOVISC contains 22 mg/mL of lightly cross-linked sodium hyaluronate (NaHA) dissolved in phosphate buffered saline. Note: The contents of the syringe are sterile; however, the product tray is non-sterile.

FOR INTRA-ARTICULAR USE. STORE AT 2°C to 25°C. PROTECT FROM FREEZING.

CAUTION: This device is restricted to sale and use by or under the supervision of a physician.

This product is for single patient use only and must not be re-sterilized. Reuse of needles or syringes used to inject this product can result in transmission of infectious agents as well as blood-borne pathogens (including HIV and hepatitis), potentially endangering patients and physicians and staff. Used needles or syringes should be discarded after each injection session and not saved for subsequent sessions on the same patient.

DO NOT USE IF PRODUCT INNER PACKAGING IS OPEN OR DAMAGED.

	European conformity and notified body number		Date of manufacture
	Sterilized using steam or dry heat		Batch, LOT number
	Single barrier system with protective packaging outside		Catalog Number
	Single Sterile barrier System		Medical Device
	Do not re-use		Do not use if package is damaged and consult instructions for use
	Do not resterilize		Importer
	Use by	 www.anikaifu.com	Consult Instructions for use
	Temperature Limit		

To report adverse events or serious incidents or opened/damaged poorly performing products, contact Anika Therapeutics, Inc. at +1-781-457-9000 or GlobalComplaints@anika.com. In case of any serious incident that has occurred in the EU, Anika Therapeutics should also contact the competent authority of your member state.

Legal Manufacturer
Anika Therapeutics, Inc
32 Wiggins Avenue
Bedford, MA 01730 USA

EU

EU Authorized Representative
Anika Therapeutics, S.r.L.
Corso Stati Uniti, 4/U
35127 Padova (PD) – Italy

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