Parcus Drill Bits for Suture Anchor Placement

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Parcus Drill Bits for Suture Anchor Placement (English)

1. Indications:

Parcus Drill Bits are provided in various lengths and diameters as complements to suture anchors and/or Drill Guides. They are used to create sockets and tunnels in bone.

2. Warnings:

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. This product must be inspected prior to each use and before sterilization.
- C. This product is not intended for use as an implant.
- D. Loading the Drill Bit in an off-axis manner may result in breakage or bending.

3. Material:

This device is manufactured from a Stainless-Steel Alloy. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. This product is made to be re-usable. However, the following inspections should be made during all stages of handling.
- B. Never use a bent Drill Bit.
- C. Never use a Drill Bit that has discoloration or staining.
- D. Never use a dull Drill Bit.
- E. Never use a Drill Bit with a faded laser mark.
- F. Never use a damaged Drill Bit.
- G. Stop use immediately if excessive force is required to drill.
- H. Never use a Drill Bit with any suspected defect or damage not covered in this IFU.

5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

6. Sterilization:

- A. This device is not provided sterile. It should be cleaned and sterilized using validated protocols before each use.
- B. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a 10^{-6} sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 Minutes	30 Minutes ¹
Pre-Vacuum	132°C / 270°F	4 Minutes	

¹ Drying times vary according to load size and should be increased for larger loads.

7. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

8. Instructions for Use:

- A. Place the Drill Bit into the chuck of a purpose designed drill or wire driver. Drill Bits may have laser etched lines or stops to help determine depth of penetration. Attach Drill Bits so that laser lines or stops are visible as required.
- B. If a Parcus Drill Guide is utilized, advance the Drill Bit through the guide and into the bone to the required depth. If the drill bit is equipped with a stop it may be used to determine the required depth.
- C. Remove the Drill Bit from the guide.
- D. Remove the Drill Guide, as necessary.

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