X-Twist Biocomposite Suture Anchor
Ancla de sutura biocompuesta X-Twist
Âncora de Sutura X-Twist Biocomposite
X-Twist Biocomposite 縫線錨釘
X-Twist مرساة الغرز المركبة أحيانيا من

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X-Twist Biocomposite Suture Anchor (EN)

1. Device Description

The X-Twist Biocomposite Suture Anchor consists of the implants (anchor and anchor tip) and the anchor driver assembly. The anchor and anchor tip are provided affixed to the driver and sterile. The X-Twist Biocomposite anchor is molded using a composite of β TCP (beta-tricalcium-phosphate) and PLGA (poly-lactic-co-glycolic acid). The anchor tip is molded using PEEK (polyetheretherketone). The anchor is fully threaded, double-helical, cannulated, and has inline fenestrations on each quarter-turn face. The anchor tip is retained on the driver via retention suture(s) or suture tape(s) that are passed through the driver cannulation, looped over the retention bridge within the tip, and returned out the proximal end of the driver handle and cleated. These devices are to be used with a drill, awl, and/or bone tap. The X-Twist Biocomposite Suture Anchors are provided sterile and in 4.75mm, 5.5mm, or 6.25mm diameters.

2. Intended Use

The X-Twist Biocomposite Suture Anchors are indicated for attachment of soft tissue to bone. These products are intended for the following indications:

Shoulder: Rotator Cuff Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoid Repair

Foot/Ankle: Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair

Knee: Anterior Cruciate Ligament Repair, Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis, Quadriceps Tendon Repair, Meniscal Root Repair, Secondary or adjunct fixation for ACL/PCL reconstruction or repair, MPFL Repair/Reconstruction

Elbow: Ulnar/Radial Collateral Ligament Reconstruction, Lateral Epicondylitis repair

3. Contraindications

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the directions for use or using the device for an indication other than those listed.

4. Intended Patient Population

Potential patients for this device include patients with symptomatic primary or recurrent tear of soft tissue amenable to repair. Patient selection factors to be considered include patient's need to obtain pain relief and improve function, patient's overall well-being, including ability and willingness to follow instructions and comply with activity restrictions.

5. Intended Users

The intended users for the device are orthopedic surgeons that are trained to perform the indicated procedures.

6. Intended Use Environments

The intended use environments for the device include hospitals, outpatient surgical settings, and ambulatory surgery sites.

7. Warnings

- A. Caution: Federal law restricts this device to sale by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Attempting to implant an anchor into hard, dense bone or into implant sites that are smaller than the recommended diameter may cause anchor failure (breakage) during insertion.
- D. Always try to approach the targeted site as close to perpendicular as possible. Off-axis anchor insertion may damage the device.

- E. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- F. The patient should be advised of the use and limitations of this device.
- G. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- H. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- I. This device must never be re-sterilized.
- J. Only the recommended instrumentation should be used to implant this device (see Section 12).
- K. This device must not be used if any of the temperature sensitive warning labels are red.

8. Cautions

- A. Applying excessive torque to the driver handle and/or paddle may cause driver failure.
- B. Improper fixation of the anchor and/or suture may leave the implant loose leading to soft tissue irritation, bone/tissue damage, and/or post operative failure.

9. Adverse Effects / Potential Complications

- A. Infection, both deep and superficial
- B. Allergies and other reactions to the device materials
- C. Loss of fixation via suture or anchor pull out can occur
- D. Suture or suture anchor breakage can occur

10. MR Safety Information

The implantable portion of these devices are MR-Safe.

11. Material Specifications

- A. Anchor Body: The anchor body is molded with a composite of β TCP (beta-tricalcium-phosphate) and PLGA (poly-lactic-co-glycolic acid).
 - B. Anchor Tip: The anchor tip is molded with PEEK (polyetheretherketone).
 - C. Shaft: Stainless Steel
 - D. Handle/paddle: Acrylonitrile Butadiene Styrene (ABS)
 - E. Suture: The suture anchors are preloaded with non-absorbable surgical sutures or suture tape made of UHMWPE (ultra-high molecular weight polyethylene). The suture or suture tape may contain undyed or dyed trace filaments of nylon or PET (polyethylene terephthalate).
- No Hazardous substances.

12. Instrumentation

Anchor Size	Bone Quality	Instrumentation
4.75mm	Soft	3.6mm Awl (XT04702)
	Normal/Hard	4mm drill (XT04705) followed by 4.75mm Tap (XT04700, XT04703)
5.5mm	Soft	4.1mm Awl (XT04701)
	Normal/Hard	4mm drill (XT04705) followed by 5.5mm Tap (XT05500, XT05501)
6.25mm	Soft	4.1mm Awl (XT04701)
	Normal/Hard	N/A – intended for soft bone only

13. Packaging and Labeling

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature (temperature indicator label is red), or has been altered in any way.
- B. The X-Twist Biocomposite Suture Anchors are provided double pouched. The device must be stored in both pouches. While the outer pouch provides a sterile and moisture barrier, the inner pouch only preserves the sterility of the device and does not prevent moisture from getting to the suture anchor. Devices stored within only one layer of packaging should be discarded.
- C. Please contact Anika Therapeutics Customer Service to report any package/label damage or alterations.

14. Sterilization

The contents were sterilized by ethylene oxide (EO).

15. Storage

These devices must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package. The X-Twist Biocomposite Suture Anchors must be stored below 25°C (77°F). Do not use the device if the temperature sensitive warning labels are red.

16. General Instructions for Use (Follow the surgical technique guide for procedure specific instructions.)

- A. Evaluate the target implant site for bone quantity and quality.
- B. Select the appropriate anchor size. The 6.25mm anchor is used primarily in bone of suboptimal quality.
- C. Inspect the device before using.
- D. Create a pilot hole using the appropriately sized awl, drill and/or tap corresponding to the implant diameter.
- E. Advance the awl into the desired anchor insertion site until laser line is flush with the bone. When preparing hard bone, advance the appropriately sized tap until the laser line is flush with the bone. Note: An improper bone preparation depth may result in poor anchor positioning.
- F. Sutures may be loaded onto the X-Twist Biocomposite Suture Anchor by passing through the distal eyelet of the anchor tip using the provided suture threader.
NOTE: The distal eyelet can accommodate up to six (6) strands of 1.6mm suture tape, #2 suture, or a combination thereof. Alternatively, the distal eyelet will accommodate up to four (4) strands of 2mm suture tape concurrently with up to three (3) #2 sutures.
- G. Align the insertion angle of the anchor tip with the prepared pilot hole and advance the anchor tip until the anchor body contacts the bone while maintaining suitable suture tension.
- H. Advance the anchor body down the driver shaft and into the prepared hole by holding the driver paddle stationary with one hand and simultaneously rotating the driver handle clockwise with the other hand until the anchor body is flush with the bone. Ensure the laser etch line on the driver is flush with the surface of the bone prior to advancing anchor body into the pilot hole.
- I. Uncleat all sutures from the handle before withdrawing the driver.
- J. Pull the driver away from the implant for removal. **Do not twist the handle upon removal.** Any remaining retention sutures may be incorporated into the repair or removed from the anchor body cannulation by pulling on one limb.
- K. Cut and remove any excess suture or suture tape limbs not used in the repair.

Reuse/Cleaning of instrumentation

Refer to the instrumentation IFU for proper cleaning, sterilization, and inspection instructions as required.

17. Device Disposal

Devices that have been removed from the packaging and not used or were explanted from a patient must be properly disposed of in accordance with the institution's policy.

18. Adverse / Reportable Event

Any serious incident that occurs in relation to the device should be reported to Anika Therapeutics at globalcomplaints@anika.com.

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